

**Identification of *Vibrionaceae* from Australian aquatic animals  
using phenotypic and PCR procedures**

J Carson MJ Higgins TK Wilson	Fish Health Unit Department of Primary Industries & Water PO Box 46, Kings Meadows Launceston, Tasmania 7249 jeremy.carson@dpiw.tas.gov.au
N Gudkovs	AAHL Australian Fish Disease Laboratory CSIRO Livestock Industries Private Mail Bag 24 Geelong Victoria 3220
TN Bryant	Medical Statistics and Computing University of Southampton Southampton General Hospital Tremona Rd, Southampton SO16 6YD Hampshire, United Kingdom

5

---

**SUMMARY**

- The *Vibrionaceae* is a large and complex group of marine bacteria that can have a significant impact on the health of aquatic animals. A range of pathogenicity is seen among the species but a consistent feature is the opportunistic basis of infection. Nearly all phases of farm production are affected from larval rearing to competent adult animals. Disease outbreaks may occur in disparate animal groups including marine mammals, fin fish, crustacea, molluscs and zooxanthellae of coral. Some strains, however, can act as probionts and have proved effective as a means of controlling disease caused by other species of *Vibrionaceae*.
- 15 Identification:** Routine, high-volume identification is achieved by phenotyping using standardised tests. To accommodate the large number of taxa and the phenotypic diversity that exists intra-species, identification is only practicable using computer-assisted probabilistic methods. The use of molecular tools for identification remains limited but PCR for several species is useful as a means of rapid screening or confirmatory identification.
- 20 Status of Australia and New Zealand:** The range of *Vibrionaceae* associated with aquatic animals in Australia is relatively small despite the diversity of habitats, geographic range and climatic variation. Major pathogens encountered are *Photobacterium damsela* ssp. *damsela*, *Vibrio anguillarum* and *Vibrio harveyi*. More unusual species isolated are *Vibrio scophthalmi* (Atlantic salmon), *Vibrio penaeicida* (southern rock lobster) and
- 25** *Photobacterium damsela* ssp. *piscicida* (southern blue fin tuna). These species are rarely encountered and appear to be incidental findings not associated with disease.

30

## Part 1. Diagnostic Overview

### Introduction

The *Vibrionaceae* is an important and ubiquitous group of bacteria in marine and estuarine environments.<sup>1</sup> Characteristically they have a close association with aquatic animals either as symbionts, normal flora or as pathogens. The ability of the group to occupy a wide diversity of habitats means that the *Vibrionaceae* have direct and indirect effects on aquaculture either as pathogens that are the cause of major disease, low grade infections that erode productivity or beneficial effects as probionts like some strains of *Vibrio alginolyticus*<sup>2</sup> or *V. mediterranei*.<sup>3</sup> To varying degrees, *Vibrionaceae* have an intimate relationship with all farmed marine species of aquatic animals from live feed inputs such as *Artemia* and rotifers, crustaceans, bivalved and single shelled molluscs to fin fish. The myriad of types, the range of species and their diverse role in aquaculture make the *Vibrionaceae* an important but challenging group of bacteria.

### Taxonomy

The *Vibrionaceae* comprise six genera: *Aliivibrio*, *Enterovibrio*, *Grimontia*, *Photobacterium*, *Salinivibrio* and *Vibrio*. The genus *Listonella* has also been proposed and although a validly published name, there is now little evidence to justify its status<sup>4</sup>, and use of the name is declining. Related to the *Vibrionaceae* is the family *Moritellaceae* and includes *M. marina* and *M. viscosa*, species pathogenic for salmonids. The taxonomy of the *Vibrionaceae* is based largely on phylogenetic analysis of small subunit ribosomal DNA, in particular, the 16S rDNA gene sequence, latterly supported by fluorescent amplified fragment length polymorphism (FAFLP) data.<sup>5,6</sup> FAFLP has proved particularly useful as a means of species delineation because of the high level of 16S rDNA sequence similarity that exists between many species of the *Vibrionaceae*. Currently some 98 species across the five genera are recognised; most of the species occur in the genus *Vibrio*. Phenotypes, when based on numerical taxonomy, are generally useful and predictive<sup>7,8</sup> but where kit systems have been used to define species, invariably descriptions are poor and insufficient for the purpose of identification. There is marked genomic plasticity within the *Vibrionaceae*, which is seen as pronounced phenotypic diversity between and within species.<sup>9,10</sup> As a consequence, finding traits that differentiate species continues to be a problem, particularly those delineated by molecular taxonomy.

### Disease in Aquatic Animals

Disease involving *Vibrionaceae* is characterised by the opportunistic basis of infection. The pathogen may have an intimate association with the host as normal flora and may exhibit its pathogenic capacity if host defences are breached. This may occur either from stress events that lead to immunosuppression, physical damage to the integument, or the emergence of aggressive biovars within a population of aquatic animals. Severity of disease may range from fulminating septicaemias, typical of disease outbreaks, through to chronic infections that affect just a few individuals.

Species of *Vibrionaceae* have an obligate requirement for NaCl<sup>11</sup> that ranges from as little as 5-15mM for *V. cholerae* and *V. metschnikovii* to 300-400 mM for *V. splendidus*<sup>12,13</sup>. A critical factor for potential pathogens is the coincidence of physiologies between host and pathogen. For osmoregulators such as fin fish, pathogens must have an optimum NaCl requirement close to physiological levels of vertebrates while for osmoconformers such as molluscs and crustacea, pathogens have optimum NaCl concentrations close to that of seawater. Extracellular pathogenicity factors such as siderophores, haemolysins, cytotoxins and proteases have been described in several species of *Vibrionaceae* and appear to be important in expression of disease.<sup>14,15,16</sup>

### Host Range

80 Development of an exhaustive list of Australian host species serves little purpose. More  
 importantly, cognisance of the range of host types is central to obtaining an understanding of  
 the pathogenic versatility of the *Vibrionaceae*. In nearly every group of marine species,  
 examples can be found of *Vibrionaceae* acting as pathogens. Major groups affected are:  
 85 marine mammals, teleost fish, crustacea, molluscs, both univalve and bivalve, and plants  
 represented by the zooxanthellae of coral.<sup>17</sup>

### Species Enzootic in Australia

Listing *Vibrionaceae* known to occur in Australia provides a guide to species that have been  
 found associated with aquatic animals, Table 1. Not all the species listed were found as  
 pathogens.

90

**Table 1.** *Vibrionaceae* enzootic in Australia

<i>Aliivibrio fischeri</i>	<i>Vibrio mediterranei</i>
<i>Photobacterium damsela</i> ssp. <i>damsela</i>	<i>Vibrio mimicus</i>
<i>Photobacterium damsela</i> ssp. <i>piscicida</i>	<i>Vibrio mytili</i>
<i>Photobacterium iliopiscarium</i>	<i>Vibrio natriegens</i>
<i>Vibrio alginolyticus</i>	<i>Vibrio navarrensis</i>
<i>Vibrio anguillarum</i>	<i>Vibrio nereis</i>
<i>Vibrio chagasii</i>	<i>Vibrio parahaemolyticus</i>
<i>Vibrio cholerae</i> non-O1	<i>Vibrio pelagius</i>
<i>Vibrio cyclitrophicus</i>	<i>Vibrio penaeicida</i>
<i>Vibrio diazotrophicus</i>	<i>Vibrio proteolyticus</i>
<i>Vibrio furnissii</i>	<i>Vibrio scophthalmi</i>
<i>Vibrio fluvialis</i>	<i>Vibrio splendidus</i> biovar I
<i>Vibrio haliotocoli</i>	<i>Vibrio tasmaniensis</i>
<i>Vibrio harveyi</i>	<i>Vibrio tubiashii</i>
<i>Vibrio ichthyenteri</i>	<i>Vibrio vulnificus</i> biovar I
<i>Vibrio lentus</i>	

### Species Exotic to Australia

Several pathogens, as listed in Table 2, have not been detected in association with aquatic  
 animals and are considered exotic.

5

**Table 2.** *Vibrionaceae* exotic to Australia

Species	Primary host
<i>Aliivibrio salmonicida</i>	salmonids
<i>Aliivibrio wodanis</i>	salmonids
<i>Moritella viscosa</i>	salmonids
<i>Vibrio ordalii</i>	salmonids
<i>Vibrio pectenocida</i>	scallops
<i>Vibrio tapetis</i>	clams
<i>Vibrio vulnificus</i> biovar II	eels

15

### Zoonotic Agents

Some species of *Vibrionaceae* are the cause of zoonoses<sup>18</sup>, Table 3. Infection invariably is  
 the result of physical trauma arising from puncture wounds or the result of ingesting  
 uncooked seafood. Life-threatening conditions have been reported occasionally in the  
 20 immunocompromised.

25

**Table 3. Zoonotic *Vibrionaceae***


---

	<i>Photobacterium damsela</i> ssp. <i>damsela</i>
	<i>Vibrio alginolyticus</i>
	<i>Vibrio cholerae</i> non-O1
	<i>Vibrio cincinnatiensis</i>
30	<i>Vibrio fluvialis</i>
	<i>Vibrio furnissii</i>
	<i>Vibrio harveyi</i>
	<i>Vibrio mimicus</i>
	<i>Vibrio parahaemolyticus</i>
35	<i>Vibrio vulnificus</i>

---

**Characteristics of the *Vibrionaceae***

Sodium chloride is critical for the growth of *Vibrionaceae*. Many species have an obligate requirement for sodium ions and growth of nearly all species is stimulated by NaCl, even those that have a low requirement for NaCl, such as *V cholerae*.<sup>19</sup> Species of the *Vibrionaceae* do not in general have fastidious growth requirements and can be readily grown on peptone-based media as long as NaCl requirements are met. Growth of *Aliivibrio salmonicida*, *A wodanis* and *M viscosa*, is more reliable on media enriched with blood. In defined media some *Vibrionaceae* require supplementation with vitamins<sup>20</sup>, while most strains, even in complex media, respond well to the addition of low levels of yeast extract.<sup>21</sup>

45 As the natural habitat of the *Vibrionaceae* is the marine environment, better growth is obtained at a slightly alkaline pH in the range of 7.5-7.8. Many species and strains will form distinctive curved rods but this characteristic is not diagnostic of the *Vibrionaceae*. In tissue smears, rods can appear preternaturally large or pleomorphic but on culture will assume more typical form and proportions.

50 Uniformly the *Vibrionaceae* are facultative anaerobes that ferment glucose. All species are oxidase positive with the exception of *V metschnikovii*, which is oxidase negative. Species are sensitive to the vibriostat 0/129, a pteridine derivative related to trimethoprim.<sup>22</sup> Some species however, such as *V lentus*, may appear resistant if inappropriate test media are used or if strains have acquired resistance to trimethoprim from the *drfA1* gene in plasmid class I integrons.<sup>23</sup> Most species of *Vibrionaceae* will grow at 25°C except for *A salmonicida*, *A wodanis*, *M marina*, *M viscosa*, *Ph iliopiscarium* and *Ph phosphoreum*, which grow at 15°C. All zoonotic species will grow at 35-37°C.

55

### Isolation Strategies

For osmoregulators, cultures from internal sites should be made on media enriched with blood such as Blood Agar Base No. 2 (Oxoid) with 7% defibrinated sheep's blood as a non-selective medium, and thiosulphate citrate bile sucrose (TCBS; Oxoid<sup>24</sup>) agar as a selective indicator medium for *Vibrionaceae*. TCBS should be used only in conjunction with a non-selective medium as not all species will grow on TCBS. For external sites on osmoregulators or any site of an osmoconformer, samples should be plated on either ZoBell's 2216E<sup>25</sup> or Johnson's marine agar<sup>26</sup> (see Appendix 1) as non-selective media and TCBS as the selective medium. For salmonids, particularly during periods of low water temperatures, samples from external sites should be plated on blood agar supplemented with 1.5% NaCl.<sup>27</sup> Plates should be incubated at 25°C for 2-3 days or at 15°C for up to 7 days for psychrophiles. Procedure for sampling should follow the ANZSDP guidelines for sample collection from finfish.<sup>28</sup>

### Preservation

Cultures of *Vibrionaceae* can be held frozen at -80°C or in liquid nitrogen. A cryopreservative (see Appendix 1) based on peptone and glycerol is suitable.<sup>29</sup> Long-term storage based on freeze-drying is effective but the menstruum must be based on meso-inositol<sup>30</sup> so as to regulate membrane phase transition temperature effects and protein stabilisation<sup>31</sup> during freeze-drying and rehydration. Some species, notably *A salmonicida* and *M marina* have proved refractory to freeze-drying and are best preserved frozen. Generally, cultures recover well from preservation. ZMA or JMA are suitable recovery media but for more fastidious species particularly *A salmonicida*, *M marina* and *M viscosa* blood agar supplemented with 2% NaCl should be used. A prudent strategy for recovery from freeze-drying is to use *Vibrio* Recovery Medium (Appendix 1), which contains sodium pyruvate<sup>32,33</sup>, which has been found useful in repair of damaged cells.

### Identification Strategies

**Phenotyping:** Typically the *Vibrionaceae* exhibit a wide diversity of phenotype both between and within species. This heterogeneity in phenotype means that identification using a small number of tests either with keys or tables is unreliable. A more dependable strategy is to use a simultaneous polythetic approach<sup>34</sup> combined with computer-assisted probabilistic identification.<sup>35</sup>

An identification matrix, VibEx7 (Appendix 3), has been developed for species of *Vibrionaceae* associated with a diversity of aquatic animals in Australia from both temperate and tropical regions<sup>8</sup>. The matrix can identify 61 species and biovars and a further 25 as yet un-named protospecies of *Vibrionaceae* from aquatic animals. The panel of tests consists of 39 biochemical and 5 antibiotic sensitivity tests; details of the tests and formulations are given in Appendix 2. The biochemical tests are in conventional tube or plate format. Alternatively, the panel of tests is available commercially in miniaturised format as MicroSys<sup>®</sup> V36 (DPIW, Launceston, Tasmania).

Probabilistic identification is undertaken using an implementation of Bayes Theorem.<sup>35</sup> An identification is reached if the Willcox probability value  $P \geq 0.99$  and the modal likelihood score  $\geq 0.001$ . The Willcox probability is a measure of the most likely identification, while the modal likelihood is a measure of the goodness-of-fit of the unknown to the nominated species.<sup>36</sup> Probabilistic identification is undertaken using PIBWin<sup>37</sup>, an intuitive software package, available freely from the University of Southampton, UK at: <http://www.som.soton.ac.uk/staff/tnb/pib.htm>

**Molecular:** Identification of *Vibrionaceae* using molecular methods is controversial. Amplified fragment length polymorphism (AFLP) and multilocus sequence typing<sup>38</sup> have proved important in establishing the taxonomic structure of the *Vibrionaceae* but these procedures are not as yet suitable for routine identification purposes. PCR amplification of

110 sequence motifs characteristic of species is a practical means of molecular identification but the scope of application is limited.<sup>39</sup> The most widely described and conserved construct is the 16S rRNA gene but for the *Vibrionaceae*<sup>40</sup> sequence divergence is only 9% which greatly limits the possibilities of identifying motifs that are species unique. Sequences other than 16S rRNA have been identified, typically virulence factors, but the constancy of these targets across strains of a species is unknown and their suitability as constructs for identification purposes is questionable.

115 Recommended PCR primer sets relevant for *Vibrionaceae* of aquatic animals in Australia are for *Ph damsela* and *V harveyi*. Based on the 16S rRNA gene both primer sets are suitable for species identification. The PCR for *Ph damsela* is in multiplex format<sup>41</sup> with one primer pair specific for *Ph damsela sensu lato* and a second primer pair for the urease gene *ureC* that is specific for *Ph damsela* ssp. *damsela*. Evidence of *Ph damsela* ssp. *piscicida* is inferred by the absence of a *ureC* amplicon. The primer pair for *V harveyi*<sup>42</sup> is compromised to some extent by known cross-reaction with some, but not all, strains of *V alginolyticus*. For the primers to be truly discriminating, positive PCR reactions must be verified using at least one of the phenotypic tests listed in Table 4.

125 **Table 4.** Differential phenotypic tests for *Vibrio harveyi* and *Vibrio alginolyticus* for confirming a positive PCR reaction for 16S rRNA *V harveyi* primers.

Species	PNPG*	Aesculin <sup>§</sup>	Putrescine <sup>□</sup>
<i>V alginolyticus</i>	16%	10%	100%
<i>V harveyi</i> bv I	95%	95%	2%
<i>V harveyi</i> bv II	100%	100%	0%

\*PNPG: 2-nitrophenyl β-D-galactopyranoside;

§Aesculin hydrolysis;

□Utilisation of putrescine

135 Identification by PCR should be limited to pure cultures and used as a means of confirming the identity of strains with atypical phenotypes. Performance of the primers has not been validated for the purpose of direct detection in tissues or environmental samples.

140 A range of primer sets for other constructs in the *Vibrionaceae* have been described, of which some have been critically evaluated. Of these, primers for the *cth* cytolysin/haemolysin gene of *V vulnificus*<sup>43</sup> and the *vah1* haemolysin gene of *V anguillarum*<sup>44</sup> appear robust. Use of these primers for identification in the absence of other species defining characteristics is not recommended since the frequency at which the targets occur intra-species is not known. The primers may have value however for screening purposes, particularly for strains where the target is known to occur or for establishing the presence of virulence factors.

### Quality Control

150 The most important factor determining success in identification is the use of standardised tests. Tests of different format should not be used unless extensive testing has been undertaken to verify test equivalence. It is important to recognise and identify sources of error that, if not well controlled, can result in unreliable identification outcomes.<sup>34</sup> Intrinsic error is associated with some tests and some species that can result in variable test outcomes. A second form of error arises from procedural deficiencies, particularly interpretation of weak positive tests. Regular use of quality control organisms (Table 5) is recommended, together with trend analysis to identify drift in performance. With practice, intra-laboratory test error of 2% is achievable.<sup>8,45,46</sup>

**Table 5.** *Vibrio* Quality Control Strains  
*Vibrio anguillarum* ATCC 19264<sup>1</sup>

160 *Vibrio fluvialis* NCTC 11327<sup>T</sup>  
*Vibrio mediterranei* CIP 103203<sup>T</sup>  
*Vibrio parahaemolyticus* ATCC 17802<sup>T</sup>  
*Vibrio tubiashii* NCIMB 1340<sup>T</sup>

### Limitations

165 The identification matrix VibEx7 reflects the diversity of *Vibrionaceae* associated with aquatic animals in Australia. It does not however include all the validly published species of *Vibrionaceae* as listed in Table 6 and these species will not be identified using the matrix.

Data for new species as they occur can be readily added to the VibEx7 matrix. The ability to differentiate the new species from those already in the matrix can be assessed using the IDSC tool in the PIBWin software.

170

**Table 6.** Species not in the VibEx7 matrix

175	Species	Host	Species	Host
	<i>Enterovibrio corallii</i>	Plant	<i>Vibrio diabolicus</i>	Environment
	<i>Enterovibrio norvegicus</i>	Animal	195 <i>Vibrio ezurae</i>	Animal
	<i>Moritella abyssi</i>	Environment	<i>Vibrio fortis</i>	Animal
	<i>Moritella japonica</i>	Environment	<i>Vibrio gallicus</i>	Animal
180	<i>Moritella profunda</i>	Environment	<i>Vibrio gigantis</i>	Animal
	<i>Moritella yanosii</i>	Environment	<i>Vibrio hepatarius</i>	Animal
	<i>Photobacterium aplysiae</i>	Animal	200 <i>Vibrio hispanicus</i>	Animal
	<i>Photobacterium frigidophilum</i>	Environment	<i>Vibrio kanaloae</i>	Animal
	<i>Photobacterium ganghwense</i>	Environment	<i>Vibrio neonatus</i>	Animal
185	<i>Photobacterium halotolerans</i>	Environment	<i>Vibrio neptunius</i>	Animal
	<i>Photobacterium indicum</i>	Environment	<i>Vibrio pacinii</i>	Environment
	<i>Photobacterium lipolyticum</i>	Environment	205 <i>Vibrio pomeroyi</i>	Animal
	<i>Photobacterium profundum</i>	Environment	<i>Vibrio ponticus</i>	Animal
	<i>Photobacterium rosenbergii</i>	Plant	<i>Vibrio rotiferianus</i>	Animal
190	<i>Vibrio aerogenes</i>	Environment	<i>Vibrio ruber</i>	Environment
	<i>Vibrio brasiliensis</i>	Animal	<i>Vibrio superstes</i>	Animal
	<i>Vibrio coralliilyticus</i>	Plant	210 <i>Vibrio xuii</i>	Animal
	<i>Vibrio crassostreae</i>	Animal		

## References

1. Urakawa H, Rivera ING. Aquatic Environment. In: Thompson FL, Austin B, Swings J, editors. *The Biology of the Vibrrios*. ASM Press, Washington, DC, 2006:175-189
- 215 2. Austin B, Stuckey LF, Robertson PAW, Effendi I, Griffith DRW. A probiotic strain of *Vibrio alginolyticus* effective in reducing diseases caused by *Aeromonas salmonicida*, *Vibrio anguillarum* and *Vibrio ordalii*. *J Fish Dis* 1995;18:93-96
3. Huys L, Dhert P, Robles R *et al*. Search for beneficial bacterial strains for turbot (*Scophthalmus maximus* L) larviculture. *Aquaculture* 2001;193:25-37
- 220 4. Austin B, Austin DA, Blanch AR *et al*. A comparison of methods for the typing of fish-pathogenic *Vibrio* spp. *Syst Appl Microbiol* 1997;20:89-101
5. Thompson FL, Swings J. Taxonomy of the Vibrrios. In: Thompson FL, Austin B, Swings J, editors. *The Biology of the Vibrrios*. ASM Press, Washington, DC, 2006:29-43
- 225 6. Thompson FL, Vandemeulebroecke K, Swings J. Genomic diversity amongst *Vibrio* isolates from different sources determined by fluorescent amplified fragment length polymorphism. *Syst Appl Microbiol* 2001;24:520-538
7. Bryant TN, Lee JV, West PA Colwell RR. Numerical classification of species of *Vibrio* and related genera. *J Appl Bacteriol* 1986;61:437-467
- 230 8. Carson J, Higgins MJ, Wilson T, Gudkovs N, Bryant TN. *Vibrrios of Aquatic Animals: Development of a National Standard Diagnostic Technology*. Fisheries Research and Development Corporation, Canberra, 2006:1-347
9. Thompson FL, Iida T, Swings J. Biodiversity of Vibrrios. *Microbiol Mol Biol Rev* 2004;68:403-431
10. Vandenberghe J, Thompson FL, Gomez-Gil B, Swings J. Phenotypic diversity amongst *Vibrio* isolates from marine aquaculture systems. *Aquaculture* 2003;219:9-20
- 235 11. Goudie ED, Gow JA. The taxonomic significance of the growth response to Na<sup>+</sup> by strains of *Vibrio*. *Can J Microbiol* 1995;41:930-935
12. Baumann P, Baumann L, Bang SS, Woolkalis MJ. Reevaluation of the taxonomy of *Vibrio*, *Beneckeia*, and *Photobacterium*: abolition of the genus *Beneckeia*. *Curr Microbiol* 1980;4:127-132
- 240 13. Baumann P, Furniss A, Lee J. Genus I *Vibrio* Pacini 1854, 411<sup>AL</sup>. In: Krieg NR, editor. *Bergey's Manual of Systematic Bacteriology*, Vol 1. Williams & Wilkins, Baltimore, 1984:518-538
14. Allam B, Paillard C, Ford SE. Pathogenicity of *Vibrio tapetis*, the etiological agent of brown ring disease in clams. *Dis Aquat Organ* 2002;48:221-231
15. Liu P-C, Lee K-K. Cysteine protease is a major exotoxin of pathogenic luminous *Vibrio harveyi* in the tiger prawn, *Penaeus monodon*. *Lett Appl Microbiol* 1999;28:428-430
- 245 16. Zhang X-H, Austin B. Pathogenicity of *Vibrio harveyi* to salmonids. *J Fish Dis* 2000;23:93-102
17. Ben-Haim Y, Zicherman-Keren M, Rosenberg E. Temperature-regulated bleaching and lysis of the coral *Pocillopora damicornis* by the novel pathogen *Vibrio coralliilyticus*. *Appl Environ Microbiol* 2003;69:4236-4242
- 250 18. Daniels NA, Shafaie AS. A review of pathogenic *Vibrio* infections for clinicians. *Infect Med* 2000;17:665-685
19. Singleton FL, Attwell R, Jangi S, Colwell RR. Effects of temperature and salinity on *Vibrio cholerae*. *Appl Environ Microbiol* 1982;44:1047-1058
20. Baumann P, Baumann L, Mandel M. Taxonomy of marine bacteria: the genus *Beneckeia*. *J Bacteriol* 1971;107:268-294
- 255 21. Smith SK, Sutton DC, Fuerst JA, Reichelt JL. Evaluation of the genus *Listonella* and reassignment of *Listonella damsela* (Love *et al*) MacDonell and Colwell to the genus *Photobacterium* as *Photobacterium damsela* comb. nov. with an emended description. *Int J Syst Bacteriol* 1991; 41:529-534
22. Matsushita S, Kudoh Y, Ohashi M. Transferable resistance to vibriostatic agent 2,4-diamino 6,7-diisopropylpteridine (0/129) in *Vibrio cholerae*. *Microbiol Immunol* 1984;28:1159-1162

- 260 23. Thungapathra M, Sinha KK, Chaudhuri SR *et al.* Occurrence of antibiotic resistance gene cassettes *aac(6)-Ib*, *dfrA5*, *dfrA12*, and *ereA2* in Class I integrons in non-O1, non-O139 *Vibrio cholerae* strains in India. *Antimicrob Agents Chemother* 2002;46:2948-2955
24. Nicholls KM, Lee JV, Donovan TJ. An evaluation of commercial thiosulphate citrate bile salt sucrose agar (TCBS). *J Appl Bacteriol* 1976;41:265-269
- 265 25. Oppenheimer CH, ZoBell CE. The growth and viability of sixty-three species of marine bacteria as influenced by hydrostatic pressure. *J Mar Res* 1952;9:10-18
26. Johnson PT. A new medium for maintenance of marine bacteria. *J Invertebr Pathol* 1968;11:144
27. Lunder T, Evensen Ø, Holstad G, Håstein T. 'Winter ulcer' in the Atlantic salmon *Salmo salar* Pathological and bacteriological investigations and transmission experiments. *Dis Aquat Organ* 1995;23:39-49
- 270 28. Handler J. Collection and submission of samples for investigation of diseases of fin fish. *National Aquatic Animal Health Technical Working Group – Advisory Document*. pp 39
29. Ward M, Watt P. The preservation of gonococci in liquid nitrogen. *J Clin Pathol* 1971;24:122-123
- 275 30. Redway KF, Lapage SP. Effect of carbohydrates and related compounds on the long-term preservation of freeze-dried bacteria. *Cryobiology* 1974;11:73-79
31. Leslie SB, Israeli E, Lighthart B, Crowe JH, Crowe LM. Trehalose and sucrose protect both membranes and proteins in intact bacteria during drying. *Appl Environ Microbiol* 1995;61:3592-3597
32. Brewer DG, Martin SE, Ordal ZJ. Beneficial effects of catalase or pyruvate in a most-probable-number technique for the detection of *Staphylococcus aureus*. *Appl Environ Microbiol* 1977;34:797-800
- 280 33. Mizunoe Y, Wai SN, Ishikawa T, Takade A, Yoshida S-I. Resuscitation of viable but unculturable cells of *Vibrio parahaemolyticus* induced at low temperature under starvation. *FEMS Microbiol Lett* 2000;186:115-120
34. Sneath PHA. Test reproducibility in relation to identification. *Int J Syst Bacteriol* 1974;24:508-523
- 285 35. Willcox WR, Lapage SP, Bascomb S, Curtis MA. Identification of bacteria by computer: theory and programming. *J Gen Microbiol* 1973;77:317-330
36. Gyllenberg HG, Niemi TK. New approaches to automatic identification of micro-organisms. In: Pankhurst RJ, editor. *Biological Identification with Computers*. Academic Press, London, 1975:121-136
37. Bryant TN. PIBWin – software for probabilistic identification. *J Appl Microbiol* 2004;97:1326-1327
- 290 38. Thompson FL, Gevers D, Thompson CC *et al.* Phylogeny and molecular identification of vibrios on the basis of multilocus sequence analysis. *Int J Syst Evol Microbiol* 2005;71:5107-5115
39. Nishibuchi M. Molecular Identification. In: Thompson FL, Austin B, Swings J, editors. *The Biology of the Vibrios*. ASM Press, Washington, DC, 2006:44-64
40. Ruimy R, Breittmayer V, Elbaze P *et al.* Phylogenetic analysis and assessment of the genera *Vibrio*, *Photobacterium*, *Aeromonas* and *Plesiomonas* deduced from small sub-unit rRNA sequences. *Int J Syst Bacteriol* 1994;44:416-426
- 295 41. Osorio CR, Toranzo AE, Romalde JL, Barja JL. Multiplex PCR assay for ureC and 16S rRNA genes clearly discriminates between both subspecies of *Photobacterium damsela*. *Dis Aquat Organ* 2000;40:177-183
42. Oakey HJ, Levy N, Bourne DG, Cullen B, Thomas A. The use of PCR to aid in the rapid identification of *Vibrio harveyi* isolates. *J Appl Microbiol* 2003;95:1293-1303
- 300 43. Brasher CW, DePaola A, Jones DD, Bej AK. Detection of microbial pathogens in shellfish with multiplex PCR. *Curr Microbiol* 1998;37:101-107
44. Hirono I, Masuda T, Aoki T. Cloning and detection of the hemolysin gene of *Vibrio anguillarum*. *Microb Pathog* 1996;21:173-182
- 305 45. Lapage SP, Bascomb S, Willcox WR, Curtis MA. Identification of bacteria by computer: general aspects and perspectives. *J Gen Microbiol* 1973;77:273-290
46. Sneath PHA, Johnson R. The influence on numerical taxonomic similarities of errors in microbiological tests. *J Gen Microbiol* 1972;72:377-392

- 310 47. Bascombe S, Lapage SP, Curtis MA, Willcox WR. Identification of bacteria by computer: identification of reference strains. *J Gen Microbiol* 1973;77:291-315
48. Taniguchi H, Hirano H, Kubomura S, Higashi K, Mizuguchi Y. Comparison of the nucleotide sequences of the genes for the thermostable direct hemolysin and the thermolabile hemolysin from *Vibrio parahaemolyticus*. *Microb Pathog* 1986;1:425-432
- 315 49. Furniss AL, Lee JV, Donovan TJ. *The Vibrios, Monograph Series 11*, Public Health Laboratory Service. HMSO, London, 1978:1-58
50. Thornley MJ. The differentiation of *Pseudomonas* from other Gram-negative bacteria on the basis of arginine metabolism. *J Appl Bacteriol* 1960;23:37-52
51. West PA, Colwell RR. Identification and Classification of *Vibrionaceae* - An Overview. In: Colwell RR. editor. *Vibrios in the Environment*. John Wiley, New York, 1984:285-363
- 320 52. Atlas R, Parks L. Bromocresol Purple Broth. In: *Handbook of Microbiological Media*. CRC Press, Boca Raton, 1993:158
53. Hansen GH, Sørheim R. Improved method for phenotypical characterization of marine bacteria. *J Microbiol Methods* 1991;13:231-241
- 325 54. Lowe GH. The rapid detection of lactose fermentation in Paracolon organisms by the demonstration of  $\beta$ -D-galactosidase. *J Med Lab Technol* 1962;19:21-25
55. Cowan ST. In: *Manual for the identification of medical bacteria*. 2<sup>nd</sup> edn. Cambridge University Press, Cambridge, 1974:1-238
56. Møller V. Simplified tests for some amino acid decarboxylases and the arginine dihydrolase system. *Acta Pathol Microbiol Scand* 1955;36:158-172
- 330 57. Smibert RM, Krieg NR. Phenotypic characterization. In: Gerhardt P. editor. *Methods for General and Molecular Bacteriology*. American Society for Microbiology, Washington, DC, 1994:607-654
58. Maslen LGC. Routine use of liquid urea medium for identifying *Salmonella* and *Shigella* organisms. *Br Med J* 1952;2(4783):545-546
- 335 59. Pitt T, Dey D. A method for the detection of gelatinase production by bacteria. *J Appl Bacteriol* 1970;33:687-691
60. Lam J, Mutharia L. Antigen-antibody reactions. In: Gerhardt P. editor. *Methods for General Bacteriology*. American Society for Microbiology, Washington, DC, 1994:104-132
61. MacFaddin JF. *Biochemical tests for Identification of Medical Bacteria*. 2<sup>nd</sup> edn. Williams & Wilkins, Baltimore, 1980:1-912
- 340 62. Hendrickson D. Reagents and Stains. In: Lennette EH. editor. *Manual of Clinical Microbiology*. 4<sup>th</sup> edn. American Society for Microbiology, Washington DC, 1985:1093-1107
63. Barrow GI, Feltham RKA. *Cowan and Steel's Manual for the Identification of Medical Bacteria*. 3<sup>rd</sup> edn. Cambridge University Press, Cambridge, 1993:1-331

345

## Part 2 – Test Methods

## Identification by Phenotype

*Principle*

350 Identification of *Vibrionaceae* is made using the entire panel of tests listed in Table 7 using the standardised test formulations given in Appendix 2. An identification is made by matching the phenotype of an unknown against the probability data matrix VibEx7 and determining the most likely identification using the computer software package PIBWin.

**Table 7.** Panel of tests for the identification of *Vibrionaceae*

355	Arginine dihydrolase
	Lysine decarboxylase
	Ornithine decarboxylase
	Acid (fermentation)
	Arbutin
	Mannitol
360	Salicin
	Sucrose
	Gentiobiose
	Growth
	7% NaCl
365	10% NaCl
	Acetoin
	Indole
	Alkaline phosphatase, pH 8.0
	Oxidase
370	Hydrolysis
	2-nitrophenyl $\beta$ -D-galactopyranoside
	L-glutamic acid 5-(4-nitroanilide)
	4-nitrophenyl sulfate
	Aesculin
375	Agar
	Gelatin
	Starch
	Sole carbon utilisation:
	$\alpha$ -ketoglutarate
380	Acetate
	D-alanine
	Citrate
	L-citrulline
	D-galactose
385	D-gluconate
	D-glucosamine
	D-glucose
	D-glucuronate
	Glycerol
390	L-histidine
	DL-3-hydroxybutyrate
	trans-4-hydroxy-L-proline
	DL-lactate
	D-lactose
395	Propionate

**Table 7.** Panel of tests for the identification of *Vibrionaceae* (cont.)

	Putrescine
	Succinate
	Sucrose
400	Resistance:
	0/129 10 µg
	0/129 150 µg
	Ampicillin 10 µg
405	Carbenicillin 100 µg
	Novobiocin 5 µg

*Sample Requirements*

Pure cultures less than 48 hours old should be used for the inoculum. Cultures recovered from preservation by freezing or drying must be subcultured at least twice before commencing identification.

*Test Procedures*

Perform an oxidase test by a preferred method. Observe agarolytic activity as pitting of colonies on maintenance agar such as JMA or ZMA.

**Inocula:** Prepare two inocula in 3 mL volumes of 2% saline: one to a density equal to McFarland 0.5, the other to McFarland 2.

**Inoculation of media:** Inoculate the decarboxylase test media with 100 µL of McFarland 2 density cell suspension. Inoculate the remaining liquid media with 100 µL of McFarland 0.5 density cell suspension. For the acetoin and arginine dihydrolase tests, inoculate the semi-solid media with a straight wire.

Sole carbon source media are spot inoculated with 2 µL of McFarland 0.5 suspension or with a multipoint inoculator. Maximum number of inocula on a plate should not exceed 30.

Gelatin and starch plates are spot inoculated with 2 µL of McFarland 0.5 suspension or with a multipoint inoculator. Maximum number of inocula on a plate should not exceed 6, well separated inoculum points.

Tests for arginine dihydrolase and decarboxylases are overlaid with 20-25 mm of sterile liquid paraffin.

Sensitivity tests are undertaken on Mueller-Hinton agar supplemented with 2% w/v NaCl. The medium is inoculated with the McFarland 0.5 suspension.


**Incubation:** Tests are incubated at 25°C for 48 hours. Observe tests daily and record changes. Sensitivity tests are incubated for 24 hours and the diameter of the zone of inhibition measured. Psychrophilic species are tested at 15°C for 8 days; sensitivity tests for 3-4 days. Known psychrophiles are listed in Table 8.

**Table 8.** Species requiring incubation at 15°C

435	<i>Aliivibrio logei</i>
	<i>Aliivibrio wodanis</i>
	<i>Aliivibrio salmonicida</i>
	<i>Moritella marina</i>
	<i>Moritella viscosa</i>
440	<i>Photobacterium iliopiscarium</i>
	<i>Photobacterium phosphoreum</i>


*Test Interpretation*

445 **Arginine dihydrolase:** A positive reaction is a pink/red colour; yellow to orange is negative. See Figure 1.



**Figure 1.** Arginine dihydrolase. **A:** negative; **B:** positive


450 **Decarboxylase tests:** The negative control should be yellow for the test to be valid; a positive test for ornithine or lysine decarboxylase is purple, See Figure 2.



**Figure 2.** Decarboxylase test. **A:** positive ornithine; **B:** positive lysine; **C:** negative control

455

**Acid from carbohydrates:** A dirty yellow to bright yellow colour is positive; pale purple to deep purple is negative. See Figure 3. Protein deamination may occur with prolonged incubation and may cause positive tests to appear negative due to alkaline pH shifts; ignore reversions.




460

**Figure 3.** Fermentation test. **A:** positive; **B:** negative

**NaCl tolerance tests:** Any signs of growth is a positive reaction.

465


**Acetoin:** Overlay semi-solid medium with 200 $\mu$ L of  $\alpha$ -naphthol followed by 100 $\mu$ L of KOH/creatine. Leave at room temperature for up to 30 minutes. Pale pink to red layer is positive; yellow to tan is negative. See Figure 4.



470 **Figure 4.** Acetoin (Voges-Proskauer) test. **A:** positive; **B:** negative

**Indole:** In a microtitre tray well add 100  $\mu$ L of Kovács' indole reagent to an equal volume of culture. Mix the contents of the well by careful aspiration with a pipette. If a pink/red colour is visible record as positive. See Figure 5. If in doubt, remove contents with a glass Pasteur pipette and allow the phases to separate in the pipette body; record the reaction based on the top phase only.

475




**Figure 5.** Indole test. Wells 1-5 are a range of positive reactions; well 6 is negative.

480 **Gelatin hydrolysis:** Flood plate with saturated ammonium sulphate. Any zone of clearing around inocula is positive.

**Starch hydrolysis:** Flood plate with Gram's iodine. Zones of yellow to light tan colour around inocula is positive; black to dark blue is negative.


485 **Urease:** A pink/red colour is positive; yellow is negative. See Figure 6.



**Figure 6.** Urease test. **A:** positive; **B:** negative


**Chromogens:**

- 490 *Indoxyl alkaline phosphatase*: Any blue/black colour is positive; no colour is negative. See Figure 7.



**Figure 7.** Indoxyl phosphate for alkaline phosphatase test, pH 8.0. **A:** negative; **B:** positive.


- 495 *Nitrophenol & nitroaniline chromogens*: Any yellow colour is positive; no colour is negative. See Figure 8.



**Figure 8.** Nitrophenol and nitroanilide chromogen test. Substrate: 2-nitrophenyl  $\beta$ -D-galactopyranoside (PNPG). **A:** positive; **B:** negative.

500


**Aesculin**: A brown to black colour is positive for hydrolysis; a very light tan to colourless is negative. See Figure 9.



**Figure 9.** Aesculin test. **A:** negative; **B:** positive

505

**Sole carbon source tests**: Examine the negative control plate; some strains and some species may show very slight growth even when using purified agar. See Figure 10.



510 **Figure 10.** Sole carbon source test. **A:** negative control; **B:** glucose. Test organism *Vibrio fluvialis* NCTC 11327<sup>T</sup>. Differentiation between positive and negative tests.

515 With reference to the control plate, examine growth on the remaining test plates. Growth in excess of the negative control is read as positive. If there is no growth with glucose the strain may be nutritionally fastidious. Retest with sole carbon source media containing 0.015g/L Casamino acids.

**Antibiotic sensitivity test:** The diameter of the zone of inhibition is interpreted using the data in Table 9; disc size is 6mm.

**Table 9.** Zone size interpretation for diagnostic antibiotics

Test	Resistant	Sensitive
0/129 10 µg	≤15 mm	≥16 mm
0/129 150 µg	No zone	Any zone
Ampicillin 10 µg	≤13 mm	≥14 mm
Novobiocin 5 µg	≤16 mm	≥17 mm
Carbenicillin 100 µg	≤22 mm	≥23 mm

#### *Probabilistic identification*

An identification is obtained using the software PIBWin and the probability matrix VibEx7. An identification is accepted if the Willcox probability score  $P$  is  $\geq 0.99$  and the modal likelihood score (MLS) is  $\geq 0.001$ .<sup>8</sup>

530 If  $P$  is  $\geq 0.99$  but the MLS is  $< 0.001$ , the unknown may represent an outlier of the nominated species. Differences between expected and observed reactions should be inspected to determine whether the suggested identification can be accepted. If test variance for unexpected results is between 30-70% then the identification can be accepted but should be reported as a poor fit that represents a species outlier. If observed and expected reactions for  
535 two or more tests are in complete disagreement then the suggested identification is erroneous and the unknown should be reported as unidentifiable.

Each species in the data matrix is preceded by the number of strains used to establish the phenotypic range. For taxa based on only a few strains an acceptable match will occur only if the unknown is very similar to the species description. Since the estimate of diversity is narrow, it is likely that strains of the species may appear as outliers.

540 In some instances, an unknown may appear to span two near related taxa such as biovars of a species. If the first and second most likely identifications are for related taxa then the probability values can be summed<sup>47</sup>, and if  $P_{\text{sum}}$  is  $\geq 0.99$ , the unknown can be assigned to the level of taxon complex or group.

545 Care must be exercised when accepting a computer-assisted identification. Ultimately professional judgement and microbiological sense must be used for interpretation.

## Identification by PCR

### Introduction

550 The primer sets have been evaluated and optimised for the purpose of culture identification. The protocols represent optimal conditions but minor refinements may be required to account for variation in the characteristics of different thermal cyclers.

### Extraction of DNA

*Reference DNA:* Extract and purify DNA from control strains using a QIAamp DNA mini-kit (cat. no. 51304, QIAGEN). Purified control DNA is used at a concentration of 50pg  $\mu\text{L}^{-1}$ .

555 *Sample DNA:* To 100  $\mu\text{L}$  of PCR grade 18M $\Omega$  water in a 1.5 mL microfuge tube, suspend sufficient cells to a density equivalent to McFarland 1. Hold the tube at 100°C in a dry-heat block for 10 minutes and then cool rapidly in a cool block for 5 minutes. Pellet the cells at 10,000 rcf for 5 minutes and collect the supernatant containing liberated DNA. The extracted DNA is suitable for amplification without purification.

### 560 PCR reaction volume

All PCR reactions are as 25  $\mu\text{L}$  volumes in 200  $\mu\text{L}$  thin-walled tubes.

### Standard PCR reagents

Standard reagents for the PCR primer sets are listed in Table 10.

	<b>Table 10.</b> Standard PCR reagents
565	PCR grade water, 18M $\Omega$
	50mM magnesium chloride
	10x Invitrogen Platinum <sup>®</sup> Taq buffer
	16mM dNTP stock (4mM each dNTP)
	Invitrogen Platinum <sup>®</sup> Taq DNA polymerase
570	Primers (20 $\mu\text{M}$ stock)

### Electrophoresis of amplicon

575 Amplicon is visualised by electrophoresis using 2% agarose gel containing 0.5  $\mu\text{g mL}^{-1}$  ethidium bromide and 1xTBE buffer, Table 11. Use a 100bp ladder as a comparative index of amplicon size. Gels should be run at 7 volts  $\text{cm}^{-1}$  constant voltage.

	<b>Table 11.</b> 10x Tris-Boric-EDTA buffer	
	Tris (base)	108.0 g
	Boric acid	55.0 g
	EDTA	8.3 g
580	pH	8.0

### *Photobacterium damsela*<sup>41</sup>

#### Targets

16S rRNA gene for species identification and *ureC* gene for subspecies identification.

#### Primers

Primer	Sequence 5' → 3'
Car1	gcttgaagagattcgagt
Car2	cacctcgcggtcttctgctg
Ure-5'	tccggaataggtaaagcggg
Ure-3'	cttgaatatccatctcatctgc

### 585 DNA controls

*Ph damsela* ssp *damsela* NCIMB 2184<sup>T</sup>

*Ph damsela* ssp *piscicida* NCIMB 2058<sup>T</sup>

*Ph iliopiscarium* ATCC 51760<sup>T</sup>

## Master mix

Component	Volume
Water	11.15 µL
50mM MgCl <sub>2</sub>	1 µL
10x reaction buffer	2.5 µL
16mM dNTPs	1.25 µL
20µM Car1	2 µL
20µM Car2	2 µL
20µM Ure-5'	2 µL
20µM Ure-3'	2 µL
5U µL <sup>-1</sup> Taq polymerase	0.1 µL
Template DNA	1 µL
Volume	25 µL

590

## Cycle parameters

Cycle parameters
95°C 4min.
30 cycles
95°C 1min.
60°C 1min.
72°C 40sec.
72°C 5min final extension

## Interpretation

595 *Ph damsela* ssp *damsela*. Species specific amplicon at 267bp for 16S rDNA with an additional amplicon at 448bp for the *ureC* gene that is diagnostic for the subspecies *damsela*.

*Ph damsela* ssp *piscicida*. A single amplicon at 267bp for 16S rDNA. Absence of an amplicon for *ureC* is diagnostic for the subspecies *piscicida*.

## Limitations

600 A single band of 267bp is diagnostic of *Ph. damsela* sensu lato. An identification of *Ph damsela* ssp *piscicida* is inferred if only this band is present but some caution needs to be used where an identification is reached on the basis of absence. Corroborating phenotypic evidence and complete 16S rRNA gene sequence should be obtained if the identification represents a new finding for a region or host not previously associated with *Ph damsela* ssp *piscicida*.

605

***Vibrio harveyi***<sup>42</sup>

## Target

16S rRNA gene

## Primers

Primer	Sequence 5' → 3'
VH-1	AACgAgTTATCTgAACCTTC
VH-2	gCAGCTATTA ACTACTACC

610 DNA controls

*V harveyi* ATCC 14126<sup>T</sup>*V alginolyticus* ATCC 17749<sup>T</sup>

## Master mix

Component	Volume
Water	18.15 µL
50mM MgCl <sub>2</sub>	1.5 µL
10x reaction buffer	2.5 µL
16mM dNTPs	1.25 µL
20µM VH-1	0.25 µL
20µM VH-2	0.25 µL
5U µL <sup>-1</sup> Taq polymerase	0.1 µL
Template DNA	1 µL
Volume	25 µL

## Cycle parameters

Cycle parameters
94°C 2min.
40 cycles
94°C 1min.
65°C 1min.
72°C 2min.
72°C 5min. final extension

## 615 Interpretation

A single band of 413bp is characteristic of *V harveyi*.

## Limitations

On the basis of testing undertaken, the primers are specific for strains of both *V harveyi* biovar I and II associated with blister disease in abalone. Some strains of *V alginolyticus* are known to cross-react because of sequence similarity with the primer regions. A positive finding must be corroborated by phenotype using the tests given in Table 4.

620

***Vibrio vulnificus***<sup>43</sup>

## Target

*cth* cytolysin/haemolysin gene.

## 625 Primers

Primer	Sequence 5' → 3'
L-CTH	ttccaacttcaaaccgaactatgac
R-CTH	gctactttctagcattttctctgc

## DNA controls

*V vulnificus* ATCC 27562<sup>T</sup>

*V parahaemolyticus* ATCC 17802<sup>T</sup>

## Master mix

Component	Volume
Water	16.375 µL
50mM MgCl <sub>2</sub>	1.25 µL
10x reaction buffer	2.5 µL
16mM dNTPs	1.25 µL
20µM L-CTH	1.25 µL
20µM R-CTH	1.25 µL
5U µL <sup>-1</sup> Taq polymerase	0.125 µL
Template DNA	1 µL
Volume	25 µL

630 *Cycle parameters*

Cycle parameters	
94°C	3min.
30 cycles	
94°C	1min.
57°C	1min.
72°C	3min.
72°C	5min. final extension

*Interpretation*

A single band of 205bp is characteristic of the cytolysin/haemolysin gene in *V vulnificus*.

*Limitations*

635 An identification of *V vulnificus* is inferred by the presence of the *cth* gene which appears to be specific for *V vulnificus*. It is noteworthy that the forward primer L-CTH is unique to *V vulnificus* while R-CTH is homologous with the thermolabile haemolysin gene of *V parahaemolyticus*.<sup>48</sup> Based on limited testing the primers appear specific for *V vulnificus*, however the frequency with which the gene occurs within the species is not known and may not be sufficiently reliable for the purpose of identification in the absence of *a priori* information.

640

*Vibrio anguillarum*<sup>44</sup>*Target*

*vah1* haemolysin gene.

645 *Primers*

Primer	Sequence 5' → 3'
VaH1-P1	accgatgccatcgctcaaga
VaH1-P2	ggatattgaccgaagagtca

*DNA controls*

*V anguillarum* ATCC 19264<sup>T</sup>

*V parahaemolyticus* ATCC 17802<sup>T</sup>

*Master mix*

Component	Volume
Water	14.4 µL
50mM MgCl <sub>2</sub>	0.75 µL
10x reaction buffer	2.5 µL
16mM dNTPs	1.25 µL
20µM VaH1-P1	2.5 µL
20µM VaH1-P2	2.5 µL
5U µl <sup>-1</sup> Taq polymerase	0.1 µL
Template DNA	1 µL
Volume	25 µL

650 *Cycle parameters*

Cycle parameters	
94°C	4min.
30 cycles	
94°C	30sec.
55°C	30sec.
72°C	60sec.
72°C	5min. final extension

*Interpretation*

A single band of 490bp is characteristic of the haemolysin gene of *V anguillarum*.

*Limitations*

655 The *vah1* gene appears to be specific to *V anguillarum* despite having a common ancestry with haemolysins from other species of *Vibrionaceae*. From limited testing however it appears that not all strains of *V anguillarum* possess the *vah1* gene and its use as a primary means of identification is limited. Primers for the *vah1* gene should be limited to determining the presence of the haemolysin gene in strains or as a supporting test.

660

**APPENDICES****Appendix 1****Maintenance Media***Johnson's Marine Agar*<sup>26</sup>

665	Peptone (Oxoid LP0037)	5.0 g
	Yeast extract	1.0 g
	Ferrous (II) sulphate (FeSO <sub>4</sub> ·7H <sub>2</sub> O)	0.2 g
	Sodium thiosulphate (Na <sub>2</sub> S <sub>2</sub> O <sub>3</sub> ·5H <sub>2</sub> O)	0.3 g
	Agar	12.0 g
	Aged seawater	900 mL
670	Distilled water	100 mL
	pH	7.5-7.6
	Autoclave at 121°C for 15 minutes	

*Sheep Blood Agar +2% NaCl*

675	Blood agar base no.2 (Oxoid CM0271)	40.0 g
	NaCl	15.0 g
	Distilled water	1000 mL
	pH	7.4±0.2
680	Autoclave at 121°C for 15 minutes and cool to 50°C; aseptically add 70 mL of sterile defibrinated sheep's blood, mix gently and pour as plates.	

*Vibrio Recovery Medium*<sup>8</sup>

685	Peptone (Oxoid LP0037)	5.0 g
	Yeast extract	1.0 g
	Ferrous sulphate (FeSO <sub>4</sub> ·7H <sub>2</sub> O)	0.2 g
	Sodium thiosulphate (Na <sub>2</sub> S <sub>2</sub> O <sub>3</sub> ·5H <sub>2</sub> O)	0.3 g
	Sodium pyruvate	1.0 g
	Bacteriological charcoal	2.0 g
	Aged seawater	900 mL
690	Distilled water	100 mL
	pH	7.5-7.6
	Autoclave at 121°C for 15 minutes	

	<i>ZoBell's Marine Agar 2216E</i> <sup>25</sup>	
695	Peptone (Oxoid LP0037)	5.0 g
	Yeast extract (Oxoid LP0021)	1.0 g
	Ferric (III) phosphate (FePO <sub>4</sub> ·4H <sub>2</sub> O)	0.1 g
	Aged seawater	900 mL
	Distilled water	100 mL
700	Agar	12.0 g
	pH	7.5-7.6
	Autoclave at 121°C for 15 minutes	
	<i>Cryopreservative (freezing)</i> <sup>29</sup>	
705	Proteose peptone no.3 (Difco)	1.0 g
	Glycerol	8 mL
	Distilled water	92 mL
	Dispense as 5 mL volumes and autoclave at 121°C for 15 minutes.	
	<i>Cryopreservative (freeze-drying)</i> <sup>49</sup>	
710	Nutrient Broth No.2 (Oxoid CM067B)	16 mL
	<i>meso</i> -inositol	2.4 g
	Filter sterilise and add aseptically to 32 mL of sterile inactivated horse serum. Mix and divide into 3 mL volumes. Store at -20°C.	
715		

## Appendix 2

### Culture Media for Identification

720 **Note:** Media formulations are based on the cited source references but have been modified in some instances for use with the *Vibrionaceae*. In most cases the modification relates to the addition of NaCl and supplementation with yeast extract.

The tests as formulated should not be substituted with variants. The media as described are standardised for the identification of *Vibrionaceae* using the identification matrix VibEx7.

	<i>Arginine Dihydrolase</i> <sup>50,51</sup>	
725	Peptone (Oxoid LP0037)	0.1 g
	Yeast extract (Oxoid LP0021)	0.1 g
	NaCl	1.5 g
	K <sub>2</sub> HPO <sub>4</sub>	0.03 g
	Phenol red (1% <i>aq</i> )	0.1 mL
730	L-arginine hydrochloride	1.0 g
	Agar	0.3 g
	Distilled water	100 mL
	pH	6.8
735	Dispense as 4mL volumes in 12x90 mm culture tubes. Autoclave at 121°C for 15 minutes.	

*Fermentation Test Medium*<sup>52</sup>

	Peptone (Oxoid LP0037)	10.0 g
	Yeast extract (Oxoid LP0021)	1.0 g
	Lab Lemco (Oxoid LP0029)	3.0 g
740	NaCl	15.0 g
	Bromocresol purple 0.5%*	8 mL
	Distilled water	900 mL
	pH	7.2

\*0.5 g bromocresol purple in 100 mL 50:50 v/v ethanol/distilled water

745 Dissolve the ingredients, check and adjust pH and dispense as 90 mL volumes in screw cap bottles. Sterilise by autoclaving at 121°C for 15 minutes.

Prepare filter sterilised stocks of the following sugars as 10% stocks in distilled water: arbutin, mannitol, salicin, sucrose and  $\beta$ -gentiobiose.

750 To 90 mL of sterile base add 10 mL of 10% sugar stock to give a final concentration of 1%. Dispense medium as 3 mL volumes in 12 x 90 mm diameter sterile tubes.

*Tolerance to NaCl*<sup>51</sup>

7% NaCl

	Tryptone (Oxoid L0042)	1.0 g
755	Yeast extract (Oxoid LP0021)	0.1 g
	NaCl	7.0 g
	Distilled water	100 mL
	pH	7.2

10% NaCl

	Tryptone (Oxoid L0042)	1.0 g
760	Yeast extract (Oxoid LP0021)	0.1 g
	NaCl	10.0 g
	Distilled water	100 mL
	pH	7.2

765 Autoclave at 121°C for 15 min as 100 mL volumes in sealed bottles to prevent evaporation. Dispense as 3mL volumes in 12x90 mm sterile tubes.

*Amylase*<sup>51</sup>

	Nutrient broth No. 2 (Oxoid CM0067)	2.5 g
770	NaCl	1.0 g
	Soluble starch	0.1 g
	Agar	1.5 g
	Distilled water	100 mL
	pH	7.5

775 Dissolve all the ingredients except for the agar. Warm to assist solution of the starch if required. Add the agar. Autoclave at 121°C for 10 minutes. Pour as plates.

*Acetoin (Voges-Proskauer) Test*<sup>49</sup>

780	Tryptone (Oxoid LP0042)	0.7 g
	Soya peptone (Oxoid LP0044)	0.5 g
	Yeast extract (Oxoid LP0021)	0.1 g
	Glucose	1.0 g
	NaCl	1.5 g
	Agar	0.3 g
785	Distilled water	100 mL
	pH	7.0

Dissolve all the ingredients including the agar. Dispense as 3 mL volumes in 12x90 mm tubes. Autoclave at 115°C for 10 min. Cool as butts.

790 *Gelatin*<sup>51</sup>

	Nutrient broth No. 2 (Oxoid CM0067)	2.5 g
	NaCl	1.0 g
	Gelatin	0.5 g
	Agar	1.5 g
795	Distilled water	100 mL
	pH	7.5

Dissolve all the ingredients except for the agar. Warm to assist solution of the gelatin. Add the agar.

Autoclave at 115°C for 20 min. Pour as plates.

800

*Indole*<sup>51,53</sup>

	Tryptone (Oxoid LP0042)	1.0 g
	Yeast extract (Oxoid LP0021)	0.1 g
	NaCl	1.5 g
805	L-tryptophan	0.04 g
	Distilled water	100 mL
	pH	7.5

Dissolve ingredients. Dispense as 3 ml volumes in 12x90 mm tubes. Autoclave at 121°C for 15min.

810

*Chromogen Tests*<sup>8,54</sup>

## Nutrient Base

	Tryptone (Oxoid LP0042)	4.0 g
	Yeast extract (Oxoid LP0021)	0.4 g
815	NaCl	6.0 g
	Distilled water	300 mL
	pH	7.5

Dissolve ingredients and dispense as four volumes of 75 mL. Autoclave at 121°C for 15 min.

## 820 Phosphate buffer: 0.01M, pH 7.5

## Stock 0.1M buffer

	NaH <sub>2</sub> PO <sub>4</sub> ·2H <sub>2</sub> O	0.245 g
	Na <sub>2</sub> HPO <sub>4</sub> ·12H <sub>2</sub> O	3.022 g
	Distilled water, made up to	100 mL
825	pH	7.5

- Prepare a stock of 0.01M buffer, pH 8.0
- |  |                       |       |
|--|-----------------------|-------|
|  | 0.1M phosphate buffer | 10 mL |
|  | Distilled water       | 90 mL |
- 830 Phosphate buffer: 0.01M, pH 8  
Prepare a x10 stock of 0.1M buffer
- |  |  |         |
|--|--|---------|
|  | NaH <sub>2</sub> PO <sub>4</sub> ·2H <sub>2</sub> O  | 0.083 g |
|  | Na <sub>2</sub> HPO <sub>4</sub> ·12H <sub>2</sub> O | 3.395 g |
|  | Distilled water, made up to                          | 100 mL  |
|  | pH   | 8.0     |
- 835 Prepare a stock of 0.01M buffer, pH 8
- |  |                       |       |
|--|-----------------------|-------|
|  | 0.1M phosphate buffer | 10 mL |
|  | Distilled water       | 90 mL |
- Chromogen stocks
- IXP medium
- 840
- |  |                              |        |
|--|------------------------------|--------|
|  | 3-indoxyl phosphate          | 0.08 g |
|  | 0.01M phosphate buffer, pH 8 | 25 mL  |
- Filter to sterilise. Aseptically add to 75 mL of nutrient base. Aseptically dispense as 2 mL volumes in 12x90 mm sterile tubes.
- PNPG medium
- 845
- |  |   |        |
|--|---|--------|
|  | 4-nitrophenyl $\alpha$ -D-galactopyranoside | 0.06 g |
|  | 0.01M phosphate buffer, pH 7.5              | 25 mL  |
- Filter to sterilise. Aseptically add to 75 mL of nutrient base. Aseptically dispense as 2 mL volumes in 12x90 mm sterile tubes.
- LGN medium
- 850
- |  |                                    |        |
|--|------------------------------------|--------|
|  | L-Glutamic acid 5-(4-nitroanilide) | 0.06 g |
|  | 0.01M phosphate buffer, pH 7.5     | 25 mL  |
- Filter to sterilise. Aseptically add to 75 mL of nutrient base. Aseptically dispense as 2 mL volumes in 12x90 mm sterile tubes.
- NPS medium
- 855
- |  |                                |        |
|--|--------------------------------|--------|
|  | 4-nitrophenyl sulphate         | 0.04 g |
|  | 0.01M phosphate buffer, pH 7.5 | 25 mL  |
- Filter to sterilise. Aseptically add to 75 mL of nutrient base. Aseptically dispense as 2 mL volumes in 12x90 mm sterile tubes.
- Protect media from light.
- 860
- Aesculin Hydrolysis*<sup>55</sup>
- |     |                              |        |
|-----|------------------------------|--------|
|     | Tryptone (Oxoid LP0042)      | 1.0 g  |
|     | Yeast extract (Oxoid LP0021) | 0.1 g  |
|     | NaCl                         | 1.5 g  |
| 865 | Aesculin                     | 0.1 g  |
|     | Ferric citrate               | 0.05 g |
|     | Distilled water              | 100 mL |
|     | pH                           | 7.5    |
- 870 Dissolve ingredients and dispense as 3 mL volumes in 12x90 mm tubes. Autoclave at 115°C for 10min. Protect medium from light.

*Carbon Source Utilisation Tests*<sup>13</sup>

## Inorganic nitrogenous base

## Buffer base

875	TRIS (basic)	6.1 g
	Distilled water	500 mL
	pH	7.5

Adjust pH with concentrated HCl

## Salts solution

880	NH <sub>4</sub> Cl	1.0 g
	K <sub>2</sub> HPO <sub>4</sub> ·3H <sub>2</sub> O	0.075 g
	FeSO <sub>4</sub> ·7H <sub>2</sub> O	0.028 g
	NaCl	11.7 g
	MgSO <sub>4</sub> ·7H <sub>2</sub> O	12.3 g
885	KCl	0.75 g
	Yeast extract (Oxoid LP0021)	0.015 g
	CaCl <sub>2</sub> ·H <sub>2</sub> O	1.45 g
	Distilled water	400 mL

890 Weight out all the salts and combine except for the calcium chloride. Add the distilled water; once the salts are fully dissolved add the calcium chloride.

Combine the buffer base and salts solution.

Divide the medium into 10 x 90 mL volumes. To each volume add 1.2 g of *purified* agar (Oxoid LP0028) and autoclave at 121°C for 15 min; cool to 55°C.

## 895 Carbon sources

Prepare 2% w/v or v/v concentrations in distilled water, of the carbon substrates listed; filter to sterilise.

$\alpha$ -ketoglutarate	L-histidine
Acetate	DL-3-hydroxybutyrate
D-alanine	trans-4-hydroxy-L-proline
Citrate	DL-lactate
L-citrulline	D-lactose
D-galactose	Propionate
D-gluconate	Putrescine
D-glucosamine	Succinate
D-glucose	Sucrose
D-glucuronate	Water (control)
Glycerol	

## Complete medium

900 To a 90 mL volume of cooled molten inorganic nitrogenous base, add 10 mL of carbon source. Mix well and pour as plates.

## Supplement for nutritionally fastidious strains

To the nitrogenous base add 0.015 g/L Casamino Acids (Difco, 0230-15).

*Decarboxylase Test*<sup>51,56,57</sup>

905	Difco Decarboxylase broth base (Møller)	300 mL
	Yeast extract (Oxoid LP0021)	0.3 g
	NaCl	4.5 g
	pH	6.5

- 910 Divide the base into three volumes of 100 mL each. To one volume add 1.0 g of L-lysine and to the second volume 1.0 g of L-ornithine; the third volume is a control. Check and adjust the pH to 6.5 if required. Autoclave at 121°C for 10 min. Check the pH of the media and, if required, adjust aseptically with 1N NaOH or HCl. Aseptically dispense the media as 3 mL volumes into 12x90 mm sterile tubes.

*Urease*<sup>58</sup>

- 915 Broth base
- |     |  |        |
|-----|--|--------|
|     | Peptone (Oxoid LP0037)                       | 0.1 g  |
|     | Glucose                                      | 0.1 g  |
|     | Sodium chloride                              | 1.5 g  |
|     | Na <sub>2</sub> HPO <sub>4</sub> (anhydrous) | 0.12 g |
| 920 | KH <sub>2</sub> PO <sub>4</sub> (anhydrous)  | 0.08 g |
|     | Phenol red (0.01%)                           | 4 mL   |
|     | Distilled water                              | 95 mL  |
|     | pH   | 6.8    |

Dissolve the ingredients and autoclave as a single volume at 115°C for 20 min.

- 925 Urea stock (40% w/v)
- |  |                 |       |
|--|-----------------|-------|
|  | Urea            | 8 g   |
|  | Distilled water | 20 mL |

Dissolve the urea and filter to sterilise.

- 930 Complete medium
- To the cool sterile base add 5 mL of sterile urea stock; mix. Dispense aseptically as 2 mL volumes in 12x90 mm sterile tubes.

*Shelf Life of Media*

- 935 Agar plates have a shelf life of 4 weeks and liquid media have a shelf life of 8 weeks when stored at 2-8°C.

*Saturated Ammonium Sulphate (Gelatin test)*<sup>59,60</sup>

- |  |                   |       |
|--|-------------------|-------|
|  | Ammonium sulphate | 10 g  |
|  | Distilled water   | 10 mL |
- 940 Store at room temperature.

*Coblentz Reagents (Acetoin test)*<sup>61</sup>

- Reagent 1
- |     |            |       |
|-----|------------|-------|
| 945 | α-naphthol | 0.5 g |
|     | Ethanol    | 10 mL |
- Store refrigerated in a dark bottle.
- Reagent 2
- |     |                 |        |
|-----|-----------------|--------|
| 950 | Creatine        | 0.03 g |
|     | KOH             | 4.0 g  |
|     | Distilled water | 10 mL  |
- Store refrigerated.

*Iodine (Starch test)*<sup>62</sup>

- 955 Gram's Iodine
- |  |                  |        |
|--|------------------|--------|
|  | Iodine           | 1.0 g  |
|  | Potassium iodide | 2.0 g  |
|  | Distilled water  | 300 mL |

Dissolve the potassium iodide in 20 mL of water and then add the iodine. Once dissolved, make up to 300 mL with water.

960 Store at room temperature in a dark bottle.

*Kovács' Reagent (Indole test)*<sup>63</sup>

	<i>p</i> -dimethylaminobenzaldehyde	2 g
965	pentan-1-ol ( <i>n</i> amyl alcohol)	30 mL
	Concentrated HCl	10 mL

Dissolve the aldehyde in the alcohol by gently warming at 50-55°C. Cool and slowly add the acid. Protect from light and store at 4°C.

*Note:* *iso* amyl alcohol is not the same as *n* amyl alcohol. The *iso*- form of pentanol cannot be used for Kovács' reagent.

970

Appendix 3: VibEx7 probability matrix for the identification of *Vibrionaceae*

Data as % strains positive		<i>Aeromonas . sobria</i> HG7	<i>G. hollisae</i>	<i>M. marina</i>	<i>M. viscosa</i>	<i>A. fischeri</i> biovar I	<i>A. fischeri</i> biovar II	<i>A. logei</i>	<i>A. salmonicida</i>	<i>A. wodanis</i>	<i>Ph. angustum</i>	<i>Ph. damsela</i> ssp. <i>damsela</i> biovar I	<i>Ph. damsela</i> ssp. <i>damsela</i> biovar II	<i>Ph. damsela</i> ssp. <i>piscicida</i>	<i>Ph. itiopiscarium</i>	<i>Ph. leiognathi</i>	<i>Ph. phosphoreum</i>
Test	No. strains	7	4	1	8	5	10	2	1	1	1	29	7	5	3	4	1
Arginine dihydrolase		86	1	1	57	1	1	1	1	1	1	99	99	99	99	75	99
Acid:	Arbutin	14	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1
	Mannitol	99	1	1	1	40	87	99	99	1	99	1	1	1	33	1	1
	Salicin	14	1	1	1	40	40	1	1	1	1	1	1	1	1	1	1
	Sucrose	57	1	1	1	20	1	1	1	99	1	1	1	1	33	50	1
	Gentiobiose	14	1	1	1	60	90	1	1	1	99	1	1	1	1	1	1
Growth:	7% NaCl	1	99	99	1	1	70	1	99	1	99	86	99	1	99	50	99
	10% NaCl	1	25	99	1	1	1	1	1	1	1	1	1	1	1	1	1
Amylase		99	1	1	99	40	55	1	1	99	1	27	14	20	99	1	1
Voges Proskauer (Acetoin)		86	1	1	1	1	1	1	1	1	1	99	99	99	99	99	99
Gelatinase		71	1	1	99	20	1	1	1	1	1	24	33	20	1	25	1
Indole		86	75	1	1	20	1	1	1	99	1	3	1	1	1	1	1
IXP alkaline phosphatase		86	50	1	13	80	99	99	1	1	99	99	85	20	33	25	99
PNPG $\alpha$ -D-galactosidase		29	1	99	1	1	20	1	1	1	1	3	1	20	1	1	99
LGN $\gamma$ -glutamyl transpeptidase		71	1	99	99	1	50	99	1	1	1	48	14	20	1	75	1
NPS sulphatase		14	1	99	1	60	99	99	1	1	99	1	1	1	1	25	1
Aesculin hydrolysis		71	1	1	20	80	99	99	1	1	50	31	1	1	1	25	1
Utilisation: $\alpha$ -ketoglutarate		14	50	99	1	1	1	99	1	1	99	3	1	20	1	1	99
	Acetate	29	99	99	88	80	1	50	99	99	99	93	1	60	99	50	1
	Alanine	1	99	99	25	20	1	1	1	1	1	10	1	99	1	1	1
	Citrate	43	99	1	1	1	60	99	99	99	99	1	1	1	99	1	1
	Citrulline	1	1	1	38	1	1	1	1	1	1	1	1	1	1	1	1
	Galactose	57	99	1	1	40	99	99	99	1	99	93	28	60	99	25	99
	Gluconate	99	75	99	1	1	20	99	99	1	1	3	14	20	99	99	1
	Glucosamine	99	25	1	50	20	99	99	1	1	99	99	42	60	99	99	99
	Glucuronate	1	1	1	1	1	1	1	1	1	1	1	1	1	33	1	1
	Glycerol	99	75	99	99	40	90	99	99	99	99	99	42	60	99	75	99
	Histidine	29	50	1	88	1	1	1	1	1	1	1	1	1	1	1	1
	DL-3-hydroxybutyrate	1	1	1	1	40	1	1	99	1	1	1	1	1	1	25	1
	Hydroxyproline	1	1	1	1	1	1	1	1	1	99	1	1	1	1	1	1
	DL-lactate	14	99	99	99	1	1	1	1	1	1	79	1	1	1	99	99
	Lactose	1	1	1	1	1	1	1	1	1	1	1	1	20	1	1	99
	Propionate	14	99	1	1	1	1	1	1	1	1	1	1	1	1	1	1
	Putrescine	14	1	1	1	1	1	1	99	1	1	1	1	1	99	1	1
	Succinate	86	75	99	99	40	90	99	99	99	99	89	42	40	99	75	99
	Sucrose	57	1	1	1	60	1	1	99	99	1	1	1	1	33	50	1
Oxidase		99	99	99	99	99	99	99	99	99	99	99	99	99	99	99	1
Agarolysis		1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1
Resistance:	0/129 10 $\mu$ g	99	1	99	38	1	1	1	1	1	1	1	14	1	1	1	1
	0/129 150 $\mu$ g	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1

Data as % strains positive		<i>Aeromonas . sobria HG7</i>	<i>G. hollisae</i>	<i>M. marina</i>	<i>M. viscosa</i>	<i>A. fischeri biovar I</i>	<i>A. fischeri biovar II</i>	<i>A. logei</i>	<i>A. salmonicida</i>	<i>A. wodanis</i>	<i>Ph. angustum</i>	<i>Ph. damsela</i> ssp. <i>damsela</i> biovar I	<i>Ph. damsela</i> ssp. <i>damsela</i> biovar II	<i>Ph. damsela</i> ssp. <i>piscicida</i>	<i>Ph. iliopiscarium</i>	<i>Ph. leiognathi</i>	<i>Ph. phosphoreum</i>
Test	No. strains	7	4	1	8	5	10	2	1	1	1	29	7	5	3	4	1
Ampicillin 10 µg		57	1	1	1	60	80	99	99	1	99	65	14	20	1	75	1
Novobiocin 5 µg		86	25	1	1	1	10	1	1	1	50	27	14	20	33	50	1
Carbenicillin 100 µg		17	1	1	1	99	80	99	99	1	99	89	42	20	33	75	1
Lysine decarboxylase		99	1	1	38	20	55	99	1	1	99	62	57	1	99	1	1
Ornithine decarboxylase		29	1	1	1	20	33	99	1	1	1	1	1	1	1	1	1
Urease		14	1	1	1	99	99	1	1	1	99	99	99	1	1	25	1

Data as % strains positive		<i>V. aestuarianus</i>	<i>V. agarivorans</i>	<i>V. alginolyticus</i>	<i>V. anguillarum</i>	<i>V. calviensis</i>	<i>V. campbellii</i>	<i>V. chagasii</i>	<i>V. cholerae</i>	<i>V. cincinnatiensis</i>	<i>V. cyclitrophicus</i>	<i>V. diazotrophicus</i>	<i>V. fluvialis</i>	<i>V. furnissii</i>	<i>V. gazogenes</i>	<i>V. haliotocoli</i>	<i>V. harveyi biovar I</i>
Test	No. strains	3	7	30	60	4	4	21	12	3	7	4	9	10	3	3	62
Arginine dihydrolase		99	1	3	96	99	1	95	1	1	99	99	99	99	1	1	15
Acid:	Arbutin	1	14	3	1	67	1	1	1	66	1	99	99	1	99	1	29
	Mannitol	99	99	99	98	50	50	99	99	99	85	99	99	99	33	99	94
	Salicin	1	57	20	1	75	75	5	1	66	1	99	99	1	99	1	79
	Sucrose	99	1	99	98	75	1	43	99	66	99	99	99	99	99	67	65
	Gentiobiose	1	99	1	1	75	50	1	1	33	14	1	11	1	1	1	95
Growth:	7% NaCl	67	29	99	93	25	99	95	58	99	99	99	99	99	99	67	99
	10% NaCl	1	1	99	3	1	1	10	1	66	1	50	56	90	66	33	5
Amylase		99	43	99	96	99	75	99	75	66	99	25	99	70	99	1	99
Voges Proskauer (Acetoin)		1	14	99	94	1	1	5	58	33	1	1	11	1	99	1	2
Gelatinase		67	1	99	77	50	99	99	67	33	99	25	44	30	99	33	97
Indole		99	1	99	93	99	99	95	99	1	85	99	99	99	1	33	97
IXP alkaline phosphatase		99	43	99	88	50	99	99	92	99	85	99	99	99	50	33	99
PNPG α-D-galactosidase		99	99	16	53	99	1	76	75	99	99	99	33	44	99	99	95
LGN γ-glutamyl transpeptidase		1	29	96	15	50	99	86	99	1	71	99	99	99	1	1	94
NPS sulphatase		1	14	40	1	1	1	99	8	1	71	25	11	1	1	1	60
Aesculin hydrolysis		33	99	10	63	99	99	99	22	99	99	99	88	10	99	33	95
Utilisation: α-ketoglutarate		67	1	99	36	25	99	99	99	1	99	1	99	99	1	1	97
	Acetate	67	43	96	51	99	25	99	99	33	99	99	99	90	99	33	89
	Alanine	99	14	99	76	99	50	99	67	1	99	99	99	99	33	1	99
	Citrate	99	1	96	95	99	99	95	92	99	99	99	99	99	66	1	95
	Citrulline	1	1	20	16	99	25	62	1	1	85	1	56	10	1	1	6

Data as % strains positive		<i>V. aestuarianus</i>	<i>V. agarivorans</i>	<i>V. alginolyticus</i>	<i>V. anguillarum</i>	<i>V. caviensis</i>	<i>V. campbellii</i>	<i>V. chagasti</i>	<i>V. cholerae</i>	<i>V. cincinnatiensis</i>	<i>V. cyclitrophicus</i>	<i>V. diazotrophicus</i>	<i>V. fluvialis</i>	<i>V. furnissii</i>	<i>V. gazogenes</i>	<i>V. halitocoli</i>	<i>V. harveyi</i> biovar I
Test	No. strains	3	7	30	60	4	4	21	12	3	7	4	9	10	3	3	62
	Galactose	99	71	13	71	99	1	99	83	66	85	99	99	99	1	99	73
	Gluconate	99	14	99	96	75	1	67	99	66	99	99	99	99	33	1	99
	Glucosamine	99	86	99	95	99	1	99	99	99	99	75	99	99	33	99	84
	Glucuronate	1	1	1	1	1	1	14	25	1	1	25	78	10	99	1	61
	Glycerol	99	1	99	95	99	99	99	99	99	99	1	99	99	99	1	85
	Histidine	1	1	99	93	50	1	10	67	1	71	99	89	99	1	1	8
	DL-3-hydroxybutyrate	1	1	1	1	99	1	1	1	1	14	1	78	80	1	1	1
	Hydroxyproline	1	1	99	1	75	1	10	1	1	71	1	11	1	1	1	50
	DL-lactate	99	29	99	91	99	75	99	99	66	99	99	99	99	66	99	99
	Lactose	1	1	1	1	1	1	1	1	1	28	50	1	1	1	1	1
	Propionate	67	1	99	6	99	1	95	58	1	99	75	99	99	66	1	94
	Putrescine	1	1	99	1	50	1	5	1	1	1	99	33	90	1	1	2
	Succinate	99	29	99	93	75	99	99	75	99	99	75	99	90	99	1	97
	Sucrose	99	1	99	96	75	1	43	99	66	99	99	89	99	99	67	65
	Oxidase	99	99	99	99	99	99	99	99	99	99	99	99	99	99	99	99
	Agarolysis	1	99	1	1	1	1	1	1	33	1	1	1	1	1	1	1
	Resistance: 0/129 10 µg	1	1	70	3	1	25	19	8	66	28	25	44	60	99	1	82
	0/129 150 µg	1	1	1	1	1	1	1	1	1	1	1	1	10	1	1	1
	Ampicillin 10 µg	1	1	99	91	1	99	99	25	1	1	1	33	80	1	33	97
	Novobiocin 5 µg	67	14	90	11	1	99	19	17	50	71	75	99	99	99	33	94
	Carbenicillin 100 µg	33	1	99	90	1	99	99	33	1	1	1	11	80	1	33	95
	Lysine decarboxylase	1	1	99	1	25	50	5	99	1	28	1	1	1	1	1	98
	Ornithine decarboxylase	1	1	66	1	1	1	1	99	1	28	1	11	1	1	1	99
	Urease	1	1	1	1	1	1	10	1	1	1	1	1	1	1	1	63

Data as % strains positive		<i>V. harveyi</i> biovar II	<i>V. ichthyoenteri</i> biovar I	<i>V. ichthyoenteri</i> biovar II	<i>V. lentus</i>	<i>V. mediterranei</i>	<i>V. metschnikovii</i>	<i>V. mimicus</i>	<i>V. mytili</i>	<i>V. natrigens</i>	<i>V. navarrensis</i>	<i>V. nereis</i>	<i>V. nigripulchritudo</i>	<i>V. ordalii</i>	<i>V. orientalis</i>	<i>V. parahaemolyticus</i>	<i>V. pectenicida</i>
Test	No. strains	23	6	3	9	19	5	7	5	12	4	6	3	6	6	18	2
	Arginine dihydrolase	1	1	1	99	84	99	1	99	1	1	99	1	1	99	6	99
	Acid: Arbutin	1	1	1	1	1	1	1	25	83	75	1	1	1	1	1	1
	Mannitol	83	60	99	75	99	99	86	99	99	99	67	1	1	99	99	1
	Salicin	4	1	1	1	80	1	1	99	92	50	1	99	1	17	18	1
	Sucrose	1	66	99	13	95	99	1	99	99	99	99	1	99	99	11	1
	Gentiobiose	9	1	1	1	21	1	1	99	99	25	1	50	1	1	1	1
	Growth: 7% NaCl	99	33	67	38	74	80	86	99	99	99	99	1	1	67	99	1
	10% NaCl	4	16	1	1	5	40	14	99	75	1	67	1	1	1	89	1
	Amylase	99	1	1	38	84	99	1	99	92	99	67	67	50	99	99	99

Data as % strains positive		<i>V. harveyi</i> biovar II	<i>V. ichthyoenteri</i> biovar I	<i>V. ichthyoenteri</i> biovar II	<i>V. lentus</i>	<i>V. mediterranei</i>	<i>V. metschnikovii</i>	<i>V. mimicus</i>	<i>V. mytili</i>	<i>V. natrigens</i>	<i>V. navarrensis</i>	<i>V. neris</i>	<i>V. nigripulchritudo</i>	<i>V. ordalii</i>	<i>V. orientalis</i>	<i>V. parahaemolyticus</i>	<i>V. pectenicida</i>
Test	No. strains	23	6	3	9	19	5	7	5	12	4	6	3	6	6	18	2
Voges Proskauer (Acetoin)		75	1	1	1	1	99	1	1	1	1	17	1	1	1	1	1
Gelatinase		83	16	1	63	32	60	71	1	42	75	17	99	99	67	99	1
Indole		99	1	1	63	95	60	99	1	27	99	83	67	1	99	99	1
IXP alkaline phosphatase		99	99	99	75	89	40	86	99	75	99	67	67	40	83	99	99
PNPG $\alpha$ -D-galactosidase		99	1	1	50	99	99	43	75	75	1	1	99	1	67	17	1
LGN $\gamma$ -glutamyl transpeptidase		87	1	1	38	99	1	99	60	67	50	99	67	99	83	99	1
NPS sulphatase		61	16	33	25	79	1	33	1	58	1	1	1	1	1	56	1
Aesculin hydrolysis		99	66	1	71	95	99	1	99	99	75	1	50	1	20	11	1
Utilisation: $\alpha$ -ketoglutarate		1	1	1	50	37	1	99	20	67	99	99	99	17	1	89	1
Acetate		99	83	99	88	99	60	99	99	92	75	99	67	1	67	99	1
Alanine		99	33	1	75	99	60	99	99	99	75	99	99	67	99	99	99
Citrate		99	16	1	13	99	40	99	99	92	99	99	99	99	99	94	1
Citrulline		1	1	1	1	11	1	1	20	83	1	83	1	1	1	17	1
Galactose		74	1	1	25	99	40	86	80	83	1	33	99	1	67	89	1
Gluconate		99	99	1	25	5	99	99	99	92	99	99	33	1	99	99	1
Glucosamine		96	99	99	13	99	20	99	99	99	99	99	99	50	83	99	50
Glucuronate		13	99	1	1	58	1	86	1	33	25	1	33	1	1	17	1
Glycerol		83	1	1	38	99	99	99	99	92	99	99	99	33	83	99	99
Histidine		1	1	1	1	89	1	99	99	99	99	99	99	33	17	99	1
DL-3-hydroxybutyrate		1	1	1	1	16	1	1	1	83	1	99	99	1	50	6	1
Hydroxyproline		1	1	1	1	1	1	1	1	25	1	1	1	1	83	83	1
DL-lactate		99	66	99	63	99	40	99	99	99	99	99	99	17	99	99	99
Lactose		1	1	1	1	99	40	1	20	8	1	1	99	1	1	1	1
Propionate		78	16	1	13	99	1	50	99	92	99	99	67	1	67	89	1
Putrescine		1	1	1	13	99	1	1	1	83	1	83	1	1	83	89	1
Succinate		99	1	99	50	95	60	99	80	99	99	83	99	1	99	89	50
Sucrose		4	66	99	1	99	99	1	99	99	99	99	1	83	99	11	1
Oxidase		99	99	99	99	99	1	99	80	99	99	99	99	99	99	99	99
Agarolysis		1	1	1	1	5	1	1	1	8	1	1	1	1	1	1	1
Resistance: 0/129 10 $\mu$ g		91	16	67	99	1	1	1	99	99	25	33	1	1	1	61	1
0/129 150 $\mu$ g		1	1	1	1	1	1	1	1	1	1	1	1	1	1	6	1
Ampicillin 10 $\mu$ g		1	33	1	13	1	40	14	1	8	50	1	1	50	1	94	50
Novobiocin 5 $\mu$ g		74	66	1	1	11	40	57	60	99	75	80	1	1	17	94	1
Carbenicillin 100 $\mu$ g		1	33	1	25	42	40	43	1	8	75	1	1	17	1	94	99
Lysine decarboxylase		99	1	1	1	22	60	99	1	1	1	1	1	1	1	99	1
Ornithine decarboxylase		1	1	1	13	1	1	99	1	1	1	1	1	1	1	99	1
Urease		1	1	1	1	1	1	1	1	1	1	1	1	1	1	17	1

Data as % strains positive		<i>V. pelagius</i> biovar I	<i>V. pelagius</i> biovar II	<i>V. penaeicida</i>	<i>V. proteolyticus</i>	<i>V. rumoiensis</i>	<i>V. scophthalmi</i>	<i>V. splendidus</i> biovar I	<i>V. splendidus</i> biovar II	<i>V. tapetis</i>	<i>V. tasmaniensis</i>	<i>V. tubiashii</i>	<i>V. vulnificus</i> biovar I	<i>V. vulnificus</i> biovar II	Phenon 6	Phenon 8	Phenon 10
Test	No. strains	10	5	4	12	3	11	37	3	1	8	7	7	5	7	7	6
Arginine dihydrolase		1	1	1	99	99	9	99	99	1	99	99	1	1	99	99	17
Acid:	Arbutin	1	1	1	1	1	1	1	1	1	1	1	71	99	99	1	1
	Mannitol	99	99	1	99	33	10	97	99	1	71	99	71	1	99	99	99
	Salicin	1	1	1	1	1	1	8	1	1	25	99	57	99	99	1	33
	Sucrose	99	80	1	1	99	99	8	1	1	25	99	1	1	86	99	83
	Gentiobiose	1	40	1	1	1	1	1	1	1	1	99	28	40	40	1	1
Growth:	7% NaCl	80	80	1	99	67	36	71	33	1	75	57	85	20	99	99	99
	10% NaCl	1	1	1	99	33	1	11	1	1	1	1	1	1	99	99	83
Amylase		20	80	33	99	67	1	84	99	99	13	99	99	99	99	71	83
Voges Proskauer (Acetoin)		1	1	1	99	1	1	1	1	1	1	1	1	1	14	1	67
Gelatinase		20	80	1	99	67	1	89	99	1	13	99	99	80	14	99	99
Indole		1	99	1	75	1	1	95	66	1	63	99	99	1	99	86	83
IXP alkaline phosphatase		99	99	50	92	99	82	92	99	1	99	86	85	99	99	99	99
PNPG $\alpha$ -D-galactosidase		99	99	99	1	67	9	92	1	1	1	99	99	80	14	29	50
LGN $\gamma$ -glutamyl transpeptidase		80	99	1	75	99	18	50	99	1	13	99	1	1	99	86	67
NPS sulphatase		70	60	1	1	1	1	79	1	1	99	1	1	1	86	86	50
Aesculin hydrolysis		99	99	1	1	33	99	97	99	1	99	99	42	25	99	1	83
Utilisation: $\alpha$ -ketoglutarate		20	20	99	99	67	9	99	99	99	38	1	99	80	99	99	99
	Acetate	90	80	99	99	67	45	76	99	1	99	86	85	99	99	99	99
	Alanine	99	99	99	99	99	1	97	99	1	99	99	99	99	99	99	99
	Citrate	80	99	99	99	99	91	97	99	99	99	99	99	99	99	99	99
	Citrulline	60	99	1	25	1	1	24	1	1	1	99	1	1	71	14	17
	Galactose	99	99	99	1	99	36	92	1	1	1	99	99	99	1	43	67
	Gluconate	99	99	99	99	67	82	84	1	1	99	99	99	99	99	99	99
	Glucosamine	70	99	99	99	99	99	89	66	99	99	99	99	80	99	86	99
	Glucuronate	1	1	25	1	1	64	21	1	1	1	86	85	80	1	1	17
	Glycerol	99	99	99	99	99	1	99	99	99	99	99	85	99	99	99	99
	Histidine	80	1	99	99	33	1	3	1	1	1	99	28	1	86	99	83
	DL-3-hydroxybutyrate	20	1	25	1	1	18	1	1	1	1	99	1	1	99	1	1
	Hydroxyproline	1	1	1	99	1	1	1	1	1	1	86	1	1	1	1	99
	DL-lactate	99	99	99	99	67	99	97	66	99	99	99	99	99	99	99	99
	Lactose	10	20	99	1	33	1	21	1	1	1	99	1	1	1	1	1
	Propionate	60	80	75	99	33	18	39	33	1	13	86	99	80	99	86	99
	Putrescine	99	99	1	99	1	1	5	1	1	1	29	1	1	99	14	17
	Succinate	99	99	99	99	99	91	99	99	99	99	86	99	99	99	86	99
	Sucrose	99	80	1	8	99	99	3	33	1	25	99	1	1	86	99	83
Oxidase		99	99	99	99	99	99	99	99	99	99	99	99	99	99	99	99
Agarolysis		1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1
Resistance:	0/129 10 $\mu$ g	60	1	25	67	1	9	29	99	1	63	14	1	1	1	57	99
	0/129 150 $\mu$ g	1	1	1	1	1	1	3	1	1	1	1	1	1	1	1	1
	Ampicillin 10 $\mu$ g	10	1	1	58	1	1	3	99	1	1	1	1	1	1	57	99
	Novobiocin 5 $\mu$ g	70	80	1	99	1	27	3	33	1	1	57	57	20	99	57	99
	Carbenicillin 100 $\mu$ g	20	1	1	1	33	1	18	99	1	1	1	1	1	1	71	99
Lysine decarboxylase		1	1	1	8	33	1	1	1	1	1	1	99	60	29	1	83
Ornithine decarboxylase		1	1	1	1	1	1	1	1	1	1	1	99	1	29	1	99
Urease		1	1	1	1	33	1	1	1	1	1	1	1	1	14	1	1

Data as % strains positive		Phenon 52	Phenon 15	Phenon 19	Phenon 20	Phenon 21	Phenon 24	Phenon 25	Phenon 26	Phenon 27	Phenon 29	Phenon 36	Phenon 41	Phenon 42	Phenon 43	Phenon 45	Phenon 46
Test	No. strains	3	8	4	6	8	4	14	16	17	17	6	6	4	11	4	6
Arginine dihydrolase		33	99	99	83	88	99	99	99	99	99	50	99	99	99	99	99
Acid:																	
	Arbutin	1	99	1	83	99	1	21	1	1	1	1	1	1	9	1	1
	Mannitol	1	99	99	67	99	99	77	99	99	88	50	1	50	91	25	99
	Salicin	1	1	1	99	99	1	57	1	1	99	1	1	1	27	1	1
	Sucrose	1	99	99	99	99	99	99	99	99	12	1	99	99	82	25	1
	Gentiobiose	1	86	1	83	99	25	57	1	1	1	1	1	1	9	1	1
Growth:	7% NaCl	99	63	25	83	99	99	43	99	88	88	50	83	75	91	25	99
	10% NaCl	67	1	1	1	1	1	1	31	1	1	17	1	1	55	1	1
Amylase		67	63	99	99	99	99	99	93	94	99	67	99	99	73	50	99
Voges Proskauer (Acetoin)		67	25	1	1	1	1	1	1	1	1	1	17	99	9	1	1
Gelatinase		99	13	1	99	99	75	99	99	99	94	1	67	50	55	25	99
Indole		1	99	99	99	88	1	93	99	99	99	50	33	1	91	25	99
IXP alkaline phosphatase		67	99	25	83	99	99	99	99	99	99	1	83	99	99	99	99
PNPG $\alpha$ -D-galactosidase		33	99	99	17	99	99	99	99	88	99	99	99	1	1	1	1
LGN $\gamma$ -glutamyl transpeptidase		1	50	99	50	99	1	14	99	82	94	67	99	99	55	25	17
NPS sulphatase		1	88	1	17	88	99	93	99	99	99	1	1	1	55	50	50
Aesculin hydrolysis		99	99	99	67	99	1	99	99	99	99	40	50	1	99	1	99
Utilisation: $\alpha$ -ketoglutarate		1	38	1	99	99	25	99	99	69	99	99	1	1	9	25	99
	Acetate	67	38	75	99	88	99	57	99	88	99	99	99	50	73	99	99
	Alanine	33	13	99	83	99	50	99	99	94	99	83	99	99	99	99	1
	Citrate	33	88	99	83	99	50	99	99	99	94	99	67	75	91	99	1
	Citrulline	1	1	25	40	88	1	1	99	1	99	50	67	25	55	1	1
	Galactose	1	99	99	83	99	99	99	99	88	99	99	99	1	9	50	99
	Gluconate	33	50	1	83	99	99	79	99	94	99	33	83	25	99	75	99
	Glucosamine	67	99	99	99	99	99	99	99	94	99	99	1	99	99	99	99
	Glucuronate	1	13	1	1	1	99	1	6	53	99	17	1	1	1	1	1
	Glycerol	67	99	99	99	99	99	99	99	99	99	99	99	99	99	99	99
	Histidine	99	1	1	33	13	1	1	18	12	1	67	17	1	9	1	1
	DL-3-hydroxybutyrate	33	1	75	67	1	1	1	1	6	6	33	1	1	1	1	1
	Hydroxyproline	33	1	1	99	1	1	1	1	1	1	1	83	1	18	1	1
	DL-lactate	33	99	99	99	99	25	99	99	94	99	99	99	99	99	75	99
	Lactose	1	1	25	1	1	1	1	1	6	6	1	1	1	1	1	1
	Propionate	67	1	99	67	99	1	7	68	88	94	83	99	99	91	75	1
	Putrescine	1	99	1	17	1	1	1	6	1	1	33	1	1	9	1	1
	Succinate	1	99	99	83	99	99	93	99	99	99	99	99	99	91	99	99
	Sucrose	99	99	99	99	99	99	99	99	99	18	17	99	99	82	25	1
Oxidase		99	99	99	99	99	99	99	99	99	99	99	99	99	99	99	99
Agarolysis		1	13	1	1	1	1	1	1	1	6	1	1	1	9	1	1
Resistance:																	
	0/129 10 $\mu$ g	67	1	1	1	1	1	14	1	12	1	1	1	1	1	1	33
	0/129 150 $\mu$ g	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1
	Ampicillin 10 $\mu$ g	33	1	1	1	1	1	1	1	1	1	1	99	75	9	25	1
	Novobiocin 5 $\mu$ g	99	1	50	17	1	1	1	18	1	1	1	67	99	18	1	1
	Carbenicillin 100 $\mu$ g	1	1	1	1	13	1	21	6	1	65	1	99	99	9	25	1
Lysine decarboxylase		1	25	1	1	1	1	1	1	1	1	1	1	1	9	1	1
Ornithine decarboxylase		1	1	1	1	1	1	1	1	1	1	1	1	1	9	1	1
Urease		1	1	1	50	1	1	1	6	1	1	1	1	1	1	1	1

Data as % strains positive		Phenon 53	Phenon 57	Phenon 58	Phenon 59	Phenon 69	Phenon 83
Test	No. strains	9	4	4	14	5	6
Arginine dihydrolase		11	99	25	1	80	99
Acid:	Arbutin	1	1	1	1	1	1
	Mannitol	99	99	99	1	99	1
	Salicin	13	1	1	7	20	1
	Sucrose	99	1	1	14	99	1
	Gentiobiose	1	25	1	93	1	1
Growth:	7% NaCl	89	75	25	21	20	1
	10% NaCl	11	1	1	1	1	1
Amylase		1	99	1	1	80	99
Voges Proskauer (Acetoin)		1	1	1	1	99	1
Gelatinase		11	75	1	21	1	99
Indole		99	75	1	7	99	1
IXP alkaline phosphatase		99	99	99	93	1	33
PNPG $\alpha$ -D-galactosidase		1	99	99	99	80	99
LGN $\gamma$ -glutamyl transpeptidase		33	1	1	1	1	66
NPS sulphatase		44	99	50	99	20	16
Aesculin hydrolysis		56	50	1	93	1	1
Utilisation: $\alpha$ -ketoglutarate		1	25	1	1	80	33
	Acetate	99	50	50	7	99	99
	Alanine	1	99	1	1	80	1
	Citrate	89	75	99	99	99	50
	Citrulline	1	1	1	1	1	1
	Galactose	1	99	99	99	60	83
	Gluconate	1	75	25	7	99	1
	Glucosamine	99	99	99	99	20	99
	Glucuronate	44	75	1	1	1	1
	Glycerol	44	99	99	99	99	99
	Histidine	1	1	1	1	60	50
	DL-3-hydroxybutyrate	89	1	1	1	60	1
	Hydroxyproline	11	1	1	1	1	1
	DL-lactate	99	75	1	1	60	99
	Lactose	1	75	99	99	20	1
	Propionate	99	1	1	1	40	1
	Putrescine	1	1	1	1	99	1
	Succinate	78	99	99	93	99	99
	Sucrose	99	1	1	7	80	1
Oxidase		99	99	99	99	99	99
Agarolysis		1	1	1	1	1	1
Resistance:	0/129 10 $\mu$ g	89	1	1	1	40	1
	0/129 150 $\mu$ g	1	1	1	1	1	1
	Ampicillin 10 $\mu$ g	11	75	50	14	1	1
	Novobiocin 5 $\mu$ g	89	1	1	1	50	1
	Carbenicillin 100 $\mu$ g	11	99	99	71	1	1
Lysine decarboxylase		1	99	99	93	1	1
Ornithine decarboxylase		1	1	1	1	1	1
Urease		1	75	99	99	80	1