



BioControl Assurance EIA - AOAC 992.11

SCOPE

This method is applicable to all foods when modified for raw meat or heavily contaminated samples.

PRINCIPLES

The detection of *Salmonella* spp. in raw meat samples is broken down into stages as follows:

- **Pre-enrichment in non-selective liquid medium**
A 1:10 dilution of the sample must be pre-enriched in buffered peptone water supplemented with 4 ml of 0.1% novobiocin solution¹ at 35-37°C for 18 to 26 h. Buffered peptone water should be warmed to room temperature or to 36 °C for large volumes. For carcass sponges, buffered peptone supplemented with 0.1% novobiocin solution is added to the moistened sponge to bring the total volume to 60-100 ml and the sample incubated at 35-37°C for 18-26 h. In the case of sponges BPW need not be warmed to room temperature before being used to re-hydrate the sponge, for all subsequent additions BPW should be warmed to room temperature.
- **Selective enrichment**
Culture from the pre-enrichment broth is inoculated into Selenite Cystine (SC) broth (35-37°C) and tetrathionate (TT) broth (42 ± 0.5°C) and incubated at the indicated temperature for 18-24h.
- **Post enrichment**
Cultures from selective liquid media are combined and inoculate into M broth and incubate at 42 ± 0.5°C for 6-8 h.
- **Enzyme immunoassay**
Follow the manufacturer's instructions retaining M broth for confirmation of presumptive positive results.
- **Cultural confirmation**
Presumptive positives can be confirmed from the retained M broth by streaking onto XLD, HE and BS agar². Confirmation carried out at an 'off-site' laboratory must be from retained BPW enrichment. Typical colonies are confirmed as *Salmonella* using biochemical and serological tests as outlined in AS 5013.10.

¹ Suspend 0.1g of novobiocin sodium salt in 100 ml of purified water. Filter sterilise (0.2 µm). Solution is stable up to 60-days in a dark bottle at 2-8 °C.

² Xylose Lysine Deoxycholate (XLD), Hektoen enteric (HE), Bismuth Sulphite (BS)

CHECKLIST

| | | |
|------------------------------|--|-------|
| Pre-enrichment | Is the buffered peptone water warmed to room temperature (to $36 \pm 1^\circ\text{C}$ for large quantities)? | _____ |
| | Is novobiocin added to BPW? | _____ |
| | Is the correct amount of enrichment broth used for the weight of sample analysed? | _____ |
| | What volume is used for carcase swabs? | _____ |
| | Is a positive control run with each batch of samples analysed? | _____ |
| | Are reference cultures inoculated into primary enrichment broth at a level of 10 to 100 cells? | _____ |
| | Is pre-enrichment done at $36 \pm 1^\circ\text{C}$ for 18-26 h? | _____ |
| Selective-enrichment | Are selective enrichment broths incubated at the appropriate temperature? | _____ |
| | Are both selective broths used on the day of preparation? | _____ |
| Post-enrichment | Is M broth incubated at $42 \pm 0.5^\circ\text{C}$ for 6-7 h? | _____ |
| | Is M broth retained for confirmation of presumptive positive samples? How is it stored? | _____ |
| | Are selective enrichments cultures combined for analysis? | _____ |
| Enzyme immunoassay | Are the manufacturer's instructions available? | _____ |
| | Are reagents stored at $2-8^\circ\text{C}$? | _____ |
| | Is incubation carried out at $35-37^\circ\text{C}$? | _____ |
| | Are kit reagents or components from other batches used in the analysis of samples? | _____ |
| Cultural confirmation | Is <i>Salmonella</i> isolated in-house from M broths? | _____ |
| | Are XLD, HE and BS agars used for confirmation? | _____ |
| (if applicable) | If an external laboratory is used is it department approved? | _____ |
| | BPW should be supplied to off-site laboratories for confirmation following AS 5013.10 | _____ |
| | Are <i>Salmonella</i> confirmed using AS 5013.10 (with regard to biochemical and serological tests)? | _____ |