



Assurance GDS for *Escherichia coli* O157:H7 in Selected Foods - AOAC 2005.04 and Assurance GDS *E. coli* O157:H7 Tq

SCOPE

These methods are applicable for detection of *E. coli* O157:H7 in meat and meat products.

PRINCIPLES

GDS is a DNA based amplification assay. Initially the organisms are allowed to grow in a proprietary modified enrichment medium (mEHEC). Any O157 present are then concentrated by using a concentration device (immunomagnetic separation), followed by amplification of highly conserved DNA sequences in the target organisms with specific primers using an automated gene amplification system. Assurance GDS^{TD} *E. coli* Tq is a modified version of AOAC 2005.04 whereby the polymerase enzyme, instead of being added separately by the analyst, is lyophilized with the other PCR reagents inside the amplification tube at the time of manufacturing.

Detection of *E. coli* O157:H7 involves the follow steps:

▪ **Enrichment**

Samples are diluted 1:10 in pre-warmed ($42 \pm 0.5^\circ\text{C}$) mEHEC medium. This method is also validated for 375 g composite samples diluted in 1.2 L or 1.5 L of mEHEC medium. Incubation is carried out for 8-18 h at $42 \pm 0.5^\circ\text{C}$. A positive control culture must be run through all procedures daily or when testing is carried out.. The sample and enrichment broth must be at the enrichment temperature for at minimum of 8h.

Note: Analysis of samples for O157 and Top 6 STEC should follow the protocols detailed in the GDA MPX Top 7 STEC methodology.

▪ **Immuno-concentration and PCR Assay AOAC 2005.04**

Enriched samples (1mL) are transferred into a sample well containing 20 μL concentrating reagent agitated at 600 rpm for 5 min. Samples are transferred using the PickPen and washed in 35 μL of resuspension buffer. The washed bead-bacteria complex (20 μL) is transferred into amplification tubes containing 10 μL of polymerase buffer and loaded in an Assurance Rotor-GeneTM (follow the manufacturer's recommended protocol).

Assurance GDS *E. coli* O157:H7 Tq

The Tq format will involve the following procedural changes:

- the volume of resuspension buffer has changed from 35 μL to 45 μL ;
- the addition of 10 μL of polymerase to the Amplification Tubes has been eliminated; and
- the volume of sample transferred from the resuspension plate to the Amplification Tubes has changed from 20 μL to 30 μL .

▪ **Interpretation**

Upon completion of the assay the Rotor-Gene program will provide a test result. Each test sample will be identified as positive, indicating that the test sample is positive for *E. coli* O157, negative indicating that the test sample is negative for *E. coli* O157, or “No Amp” indicating that amplification did not occur.

A “No Amp” reading may be due to reagent or test failure or operator error. In this event the test must be repeated using the same enrichment cultures. If the result continues to show “No Amp” the equipment supplier must be contacted for technical services. In this case, the enrichment broth must be analysed using an alternate method or the product deemed positive for *E. coli* O157:H7 for disposition purposes

▪ **Confirmation of Positive results**

Positive samples must be confirmed using an approved method at an approved laboratory or the product deemed positive for *E. coli* O157:H7 for the purposes of disposition.

CHECKLIST

Enrichment	Is the sample enriched in mEHEC medium?	_____
	Is enrichment carried out at 42 ± 0.5°C for 8-18 h?	_____
	Is the correct amount of enrichment broth used?	_____
	Is a positive control run with each batch of samples/daily?	_____
	Are reference cultures inoculated into enrichment media at a level of 10-100 cells/mL?	_____
Assay	Are manufacturer’s instructions available for reference (what method is being used)?	_____
	What volume of resuspension buffer is used during sample preparation?	_____
	What volume of sample wash is added to Amplification tubes?	_____
	Is polymerase prepared and added to Amplification tubes on the day of use?	_____
	Are internal controls run with each batch of samples?	_____
	Are technicians familiar with and trained in the operation of the GDS Automated System and the Rotor-Gene Program?	_____
	Is the shelf-life of media controlled?	_____
Confirmation	Is confirmation carried out from the enrichment culture?	_____
	Is confirmation carried out using an approved method at an approved laboratory?	_____