



## Reveal (20-hour) – AOAC 2000.14

### SCOPE

This method is applicable for the detection of *E. coli* O157:H7 in environmental swabs of stainless steel, raw ground beef and raw beef cubes.

### PRINCIPLES

Specific antibodies for *E. coli* O157:H7 are bound to a colloidal gold and separately to a solid support. *E. coli* O157:H7 reacts with gold conjugated antibodies forming a complex which flows across a lateral flow membrane where it is bound by antibody immobilised on the membrane. Positive reaction is indicated by a coloured line in the sample test window. A line in the test verification window indicates that the analysis has proceeded correctly

#### ▪ **Enrichment**

Place environmental swabs of stainless steel into 225 ml of modified trypticase soya broth plus novobiocin (mTSB+N, FDA BAM). Incubate for 20 h at 35 – 37°C. Cool to room temperature. For raw beef samples add 9 x sample weight of modified *E. coli* broth plus novobiocin (mEC+n). Incubate for 20 h at 35 – 37°C. Cool to room temperature. This method has been validated for use with 375g samples enriched in one litre of broth.

#### ▪ **Immunoassay**

Carry out the Reveal assay following the manufacturer's instructions, incubating at room temperature. Ensure that broths are equilibrated to room temperature prior to commencing the assay

#### ▪ **Confirmation**

The manufacturer's instructions are not to be followed for confirmation of presumptive positives. Confirmation of *E. coli* O157 is to be carried out from enrichment broths using an Immunomagnetic separation procedure ie ISO 16654:2001, FSIS MLG 5 or FDA BAM Chapter 4A. Confirmation must be undertaken out at a department approved laboratory

## CHECKLIST

<b>Enrichment</b>	Environmental swabs	_____
	Is 225 ml of mTSB+n (or Reveal 20h media) used for enrichment of samples?	_____
	Is novobiocin added on the day of use?	_____
	Is enrichment carried out at 35 - 37°C for 20 h?	_____
	Raw meat samples	_____
	Is mEC+n (or Reveal 20h media) used for enrichment of samples?	_____
	Is novobiocin added (mEC+n) on the day of use?	_____
	What weight of sample and volume of enrichment broth are used?	_____
	Is enrichment carried out at 35 - 37°C for 20 h?	_____
	Is a positive control run with each batch of samples analysed?	_____
<b>Immunoassay</b>	Are reference cultures inoculated into primary enrichment broth at a level of 10 to 100 cells?	_____
	Are the manufacturer's instructions available and are they reproduced in the laboratory manual?	_____
	Are technicians familiar with positive and negative reactions?	_____
<b>Cultural confirmation</b>	Is incubation carried out at room temperature?	_____
	Is <i>E. coli</i> O157:H7 confirmed from initial enrichment broth at a department approved laboratory?	_____
	Is an approved method used for confirmation (ISO 16654:2001, MLG 5, FDA BAM Chapter 4A)?	_____
	Is IMS concentration included in the confirmation methodology?	_____