



BACGene *E. coli* O157:H7 Workflow – AOAC 022002

SCOPE

This method is applicable to raw ground beef, raw beef trim, and other selected foods. This method has been validated in 375 g of raw ground beef and beef trim samples.

PRINCIPLES

The BACGene *E. coli* O157:H7 workflow is a qualitative real-time PCR assay that utilizes BACGene *E. coli* O157:H7 and BACGene Mplex STEC Screen test kits for detection of *E. coli* O157:H7. Target DNA fragments when detected are amplified and then identified with a fluorescence-labelled hybridization probe.

Detection of *E. coli* O157:H7 involves the follow steps:

- **Enrichment**
Sample (375 g) is enriched in 750 mL (1:2) of pre-warmed (to $41.5 \pm 1^\circ\text{C}$) mTSB. Sample and enrichment media are placed in a stomacher bag fitted with a side filter and homogenized by hand massage or using a stomacher. Incubation is carried out at $41.5 \pm 1^\circ\text{C}$ for 10 - 24 h. It is essential that the temperature of the broth and sample is at $41.5 \pm 1^\circ\text{C}$ for a minimum of 10 h. A positive and a negative control culture must be run through all procedures daily or when testing is carried out.
- **PCR Assay**
Sample preparation for bacterial DNA extraction and BACGene *E. coli* O157:H7 PCR assay is carried out following the manufacturer's recommended protocol. Samples that return an O157 positive (H7¹ negative or positive) must be tested using the BACGene Mplex STEC Screen PCR assay for confirmation of the Shiga toxin encoding genes *stx1/stx2* and *eae* as per the manufacturer's recommended protocol.
- **Interpretation**
Samples that are positive for O157 by BACGene *E. coli* O157:H7 assay and positive for *stx1/stx2* and *eae* genes by BACGene Mplex STEC Screen PCR assay are regarded as potential positive and must undergo confirmation.

If the internal positive control is invalid in negative samples, the test must be repeated using the same enrichment cultures. If the internal positive control for the re-test sample is still invalid, the enrichment sample should be diluted 1:5 and 1:10 in nuclease-free water and re-tested in PCR. If inhibition remains in the diluted samples the equipment supplier must be contacted for advice, and the enrichment broth must be analysed using an alternate method or the sample deemed positive for *E. coli* O157:H7.
- **Confirmation of positive results**
For all positive samples and negative samples with an invalid positive control result, enriched media must be analysed using a department approved confirmatory method. Confirmation must be carried out at a Department approved laboratory.

¹ Testing for the H7 gene/antigen is not required

CHECKLIST

Enrichment	Is the enrichment media pre-warmed to $41.5 \pm 1^\circ\text{C}$ before use?	_____
	Is enrichment carried out at $41.5 \pm 1^\circ\text{C}$ and is the enrichment broth and sample at $41.5 \pm 1^\circ\text{C}$ for a minimum of 10 h?	_____
	Is the correct amount of enrichment broth used?	_____
	Is a positive and a negative control run with each batch of samples/daily?	_____
	Are reference cultures inoculated into enrichment media at a level of 10-100 cells per sample?	_____
PCR Assay	Are manufacturer's instructions followed in regards of cross-contamination prevention?	_____
	Are internal controls run with each batch of samples?	_____
	Are technicians familiar with and trained in the operation of PCR automated instruments and the associated software?	_____
	Are O157 positive samples further analysed by BACGene Mplex STEC Screen PCR assay?	_____
	Is the shelf-life of media and kits controlled?	_____
Confirmation	Is confirmation carried out from the enrichment culture (mTSB)?	_____
	Is confirmation carried out using a department approved confirmatory method at a department approved laboratory?	_____