



## Detection, Isolation and Identification of Top Seven Shiga Toxin-Producing *Escherichia coli* (STEC) from Meat and Meat Products - MLG 5C.04

### SCOPE

This method is applicable for detection and isolation of Top 7 Shiga-toxin producing *E. coli* (STEC) (*E. coli* O157, O26, O45, O103, O111, O121 and O145) in meat products.

### PRINCIPLES

MLG 5C.04 utilises GENE-UP PCR detection assays, GENE-UP® STEC - *stx* & *eae* (EH1) and GENE-UP® Pathogenic *E. coli* (PEC), in the screening followed by cultural isolation. EH1 PCR detects the *stx* and *eae* genes in a sample, while PEC PCR detects the presence of *stx* and *eae* genes in the same cell. Positive samples then undergo serogroup specific PCR analysis (GENE-UP STEC Top 6 and O157:H7 PCR). Cultural isolation of the screen positive sample requires the use of immuno-concentration using VIDAS® *E. coli* serogroups (ESPT), followed by plating on modified Rainbow Agar (mRBA). Colonies are then tested for specific O antigens using latex agglutination and PCR. Positive colonies are purified on Sheep Blood Agar (SBA) for confirmation. Colonies on SBA are tested for the presence of O antigens, then identified using the Bruker® MALDI Biotyper and a confirmation rapid screening PCR.

The detection of STEC can be broken down into the following steps:

- **Enrichment**  
Samples (325 ± 32.5 g) are diluted in 975 ± 19.5 mL mTSB. If meat pieces are overweight, prepare a second sub-sample that must be ≥ 63 and ≤ 357.5 g at 1:4 dilution. Samples and diluent are stomached and incubated static at 42 ± 1°C for 15-24 h. A positive control must be included. It is recommended that a negative control (*E. coli* ATCC 25922) and a blank are also run with each batch of samples.  
  
Note: DAFF export sample weight collected is 375 ± 37.5 g. Therefore, when using this method, a sub-sample needs to be prepared and analysed.
- **Rapid Screening PCR for *stx/eae* and O-group**  
Samples are screened for the presence of *stx* and *eae* using EH1 and PEC PCR assays following the GENE-UP User's Guide. Samples negative for *stx* and *eae* targets are considered negative for Top 7 STEC. Samples that are positive for *stx* and *eae* genes but negative for PEC PCR are considered negative. Samples that test positive on the PEC PCR will be further analysed by GENE-UP STEC Top 6 and GENE-UP O157:H7 PCR to determine if one of the Top 7 serogroups is present. Samples negative for these serogroups are considered negative for Top 7 STEC. A positive result for *stx/eae* and O group is regarded as potential positive.
- **Isolation and Identification**  
**Immuno-concentration and streaking onto mRBA**  
All potential positive samples from PCR screening are immuno-concentrated on the VIDAS® *E. coli* serogroups (ESPT) kit by using the ESP2 Protocol followed by plating onto mRBA. All plates are incubated at 35 ± 2 °C for 20-24 h.  
**Examination of mRBA**  
*E. coli* O157:H7 colonies typically have black or grey colouration. For all six non-O157, see MLG 5C Appendix 2 for colony characteristics.  
**Presumptive Rapid Screening**  
Colonies are picked from all plates and tested for agglutination with O antiserum. At least one colony of each morphological type on each plate is tested using latex and *stx/eae* PCR assay. Samples that have no growth on mRBA or colonies that are agglutination/PCR negative are reported as negative for STEC. Agglutination/PCR positive samples are regarded as presumptive positive.
- **Confirmation**  
Presumptive positive colonies are to be streaked on SBA (incubate at 35 ± 2°C for 16 – 24 h) for confirmation including a positive control of interest. Check for contamination. Perform proteomic confirmatory tests using Bruker MALDI Biotyper or validated equivalent system (or biochemical<sup>1</sup> tests, see end note next page). Perform serological identification using latex agglutination. Agglutination sample must be further analysed for *stx* and *eae* genes by PCR and serogroup PCR. A sample that confirms as an *E. coli* isolate containing *stx* and *eae* genes, and one of the target serogroups is regarded as confirmed positive STEC.

**CHECKLIST**

<b>Enrichment</b>	Is the sample enriched in mTSB?	_____
	Is enrichment carried out at 42 ± 1 °C for 15-24 h?	_____
	Is the correct amount of enrichment broth used (i.e. 975 ±19.5 mL of mTSB for 325 ± 32.5 g sample)?	_____
	Is a sub-sample processed and enriched at 1:4 dilution (i.e. one portion of meat in three portions of broth)?	_____
	Are a positive control and a blank run with each batch of samples analysed?	_____
	Are control cultures inoculated into enrichment broth at a level of 10 to 100 cells?	_____
<b>Screening</b>	Is screening for <i>stx</i> and <i>eae</i> undertaken using EH1 and PEC PCR assays?	_____
	Is analysis for serogroup specific genes carried out using GENE-UP STEC Top 6 and O157:H7 PCR assays?	_____
	Is a cocktail of Top 7 STEC cultures run with each PCR?	_____
<b>Immuno-concentration</b>	Is immuno-concentration performed using the VIDAS <i>E. coli</i> Serogroups (ESPT) kit?	_____
	Is the ESP2 Protocol used for immuno-concentration?	_____
<b>Isolation</b>	Are filtrates and dilutions from immuno-concentrate used to streak on mRBA?	_____
	Are all plates incubated at 35 ± 2 °C for 20-24 h?	_____
	Are all morphological types from all plates confirmed by latex agglutination and then PCR?	_____
<b>Confirmation</b>	Are latex and PCR positive colonies streaked onto SBA and incubated at 35 ± 2 °C for 16-24 h?	_____
	Are colonies on SBA analysed for <i>E. coli</i> using a Bruker MALDI Biotyper or a validated equivalent system?	_____
	Are a portion of colonies on SBA confirmed by latex agglutination?	_____
	Are latex positive colonies confirmed for <i>stx/eae</i> by PCR assays?	_____
	Are latex positive colonies confirmed for serogroups using serogroup specific PCR?	_____

<sup>i</sup>Biochemical Confirmation

A laboratory may elect to use biochemical confirmation methods (VITEK® 2) if a Bruker® MALDI Biotyper is unavailable or there is an interruption in the reagent supply chain. Inoculate isolate onto VITEK® 2 GN cards (if using VITEK® 2 Compact) or equivalent. An isolate is positive on VITEK® 2 if biochemically identified as *E. coli*.